



For Professional Use Only

eSens Influenza virus A H5/H7/H9 QL PCR kit

REF ES3303B

Instructions for Use

1 INTENDED USE

eSens Influenza virus A H5/H7/H9 QL PCR kit is an *in vitro* nucleic acid amplification test for typing (identification) of *Influenza virus A* subtypes H5, H7, H9 in *Influenza virus* cultures and biological material containing *Influenza virus A* RNA using real-time hybridization-fluorescence detection of amplified products.

eSens Influenza virus A H5/H7/H9 QL PCR kit can be used with suspected influenza without distinction of form and presence of manifestation.

The material for PCR analysis is the cDNA samples obtained from human biological material: nasal swabs (inferior nasal meatus), oropharyngeal swabs (posterior pharyngeal wall), sputum (or tracheal aspirates), bronchoalveolar lavage, autopsy material in which the *Influenza virus A* RNA was detected. In case of lower respiratory tract diseases (bronchitis, bronchiolitis, pneumonia) the most informative material is sputum (or tracheal aspirates) and bronchoalveolar lavage.

PCR kit should be used for analysis of cDNA samples in which the *Influenza virus A* RNA was detected within the analysis of biological material and viruses cultures.

NOTE: The results of PCR analysis are taken into account in complex diagnostics of disease.

2 PRINCIPLE OF PCR DETECTION

Typing (identification) of subtypes H5, H7, H9 by the polymerase chain reaction (PCR) is based on the amplification of the haemagglutinin gene fragments of given subtypes of *Influenza virus A* using specific primers. In the real-time PCR, the amplified product is detected with the use of three fluorescent dyes. These dyes are linked to oligonucleotide probes, which bind specifically to the amplified product during thermocycling. The real-time monitoring of fluorescence intensities during the real-time PCR allows the detection of accumulating product without re-opening the reaction tubes after the PCR run.

eSens Influenza virus A H5/H7/H9 QL PCR kit uses “hot-start”, which greatly reduces the frequency of nonspecifically primed reactions. “Hot-start” is guaranteed by the separation of nucleotides and Taq-polymerase by using chemically modified polymerase (TaqF). The chemically modified polymerase (TaqF) is activated by heating at 95 °C for 15 min.

The results of amplification are registered in the following fluorescence channels:

Table 1

Channel for fluorophore	FAM	JOE	ROX
DNA-target	<i>Influenza virus A H5</i> cDNA	<i>Influenza virus A H7</i> cDNA	<i>Influenza virus A H9</i> cDNA
Target gene	Haemagglutinin gene	Haemagglutinin gene	Haemagglutinin gene

3 CONTENT

eSens Influenza virus A H5/H7/H9 QL PCR kit (ES3303B) includes:

Reagent	Description	Volume, ml	Quantity
PCR-mix-FL <i>Influenza virus A H5, H7, H9</i>	clear liquid from colorless to light lilac colour	0.6	1 tube
PCR-buffer-B	colorless clear liquid	0.3	1 tube
Polymerase (TaqF)	colorless clear liquid	0.03	1 tube
Positive Control <i>Influenza virus A H5, H7, H9</i> (C+<i>Influenza virus A H5, H7, H9</i>)	colorless clear liquid	0.2	1 tube
TE-buffer	colorless clear liquid	0.2	1 tube

eSens Influenza virus A H5/H7/H9 QL PCR kit is intended for 55 reactions (including controls).

4 ADDITIONAL REQUIREMENTS

- Disposable powder-free gloves and a laboratory coat.
- Pipettes (adjustable).
- Sterile RNase-free pipette tips with aerosol filters (up to 200 µl).
- Tube racks.
- Vortex mixer.
- PCR box.
- Real-time instruments (for example, Rotor-Gene Q (QIAGEN, Germany), CFX 96 Touch, CFX 96 Opus (Bio-Rad, USA), QuantStudio 5 (Thermo Fisher Scientific), or equivalent).
- Disposable polypropylene PCR tubes:
 - a) 0.2-ml PCR tubes with optical transparent domed or flat caps if a plate-type instrument is used;
 - b) 0.2-ml PCR tubes with flat caps or strips of four 0.1-ml Rotor-Gene PCR tubes if a rotor-type instrument is used.
- Refrigerator with the range from 2 to 8 °C.

- Deep-freezer with the range from minus 24 to minus 16 °C.
- Reservoir for used tips.

5 GENERAL PRECAUTIONS

The user should always pay attention to the following:

- Use sterile pipette tips with aerosol filters and use a new tip for every procedure.
- Store all extracted positive material (specimens, controls and amplicons) away from all other reagents and add it to the reaction mix in a distantly separated facility.
- Thaw all components thoroughly at room temperature before starting an assay.
- When thawed, mix the components and centrifuge briefly.
- Use disposable protective gloves and laboratory cloths, and protect eyes while samples and reagents handling. Thoroughly wash hands afterwards.
- Do not eat, drink, smoke, apply cosmetics, or handle contact lenses in laboratory work areas.
- Do not use a kit after its expiration date.
- Dispose of all specimens and unused reagents in accordance with local regulations.
- Samples should be considered potentially infectious and handled in biological cabinet in compliance with appropriate biosafety practices.
- Clean and disinfect all samples or reagents spills using a disinfectant, such as 0.5 % sodium hypochlorite or another suitable disinfectant.
- Avoid inhalation of vapors, samples and reagents contact with the skin, eyes, and mucous membranes. Harmful if swallowed. If these solutions come into contact, rinse the injured area immediately with water and seek medical advice if necessary.
- Safety Data Sheets (SDS) are available on request.
- Use of this product should be limited to personnel trained in DNA amplification techniques.
- Workflow in the laboratory must be one-directional, beginning in the Extraction Area and moving to the Amplification and Detection Area. Do not return samples, equipment and reagents in the area where the previous step was performed.



Some components of this kit contain sodium azide as a preservative. Do not use metal tubing for reagent transfer.

6 SAMPLING AND HANDLING

eSens Influenza virus A H5/H7/H9 QL PCR kit is intended for analysis of the RNA extracted with RNA extraction kits from *Influenza virus* cultures and biological material containing *Influenza virus A* RNA (nasal swabs (inferior nasal meatus), oropharyngeal swabs (posterior pharyngeal wall), sputum (or tracheal aspirates), bronchoalveolar lavage and autopsy material).

Pretreatment

Influenza virus cultures testing is recommended to carry out after the prior dilution to the concentration not more than 10^5 GE/ml.

7 WORKING CONDITIONS

eSens Influenza virus A H5/H7/H9 QL PCR kit should be used at 18–25 °C.

8 PROTOCOL

8.1 RNA extraction

Any commercial nucleic acid extraction kit, if IVD-CE validated for the indicated specimen types, could be used.

Ecoli Dx, s.r.o. recommends:

- For the manual extraction

- **RIBO-prep** (K2-9-Et-100-CE)

- For the automatic extraction

- **ePure Viral Nucleic Acid Extraction Kit** (E2003)

NOTE: Extract the RNA according to the manufacturer's protocol.

8.2 Reverse transcription

It is recommended to use the following kit for the complementary DNA (cDNA) synthesis from the RNA:

- **REVERTA-L** (K3-4-50-CE); (K3-4-100-CE)

NOTE: Carry out the reverse transcription according to the manufacturer's protocols.

8.3 Preparing PCR

8.3.1 Preparing tubes for PCR

The type of tubes depends on the PCR instrument used for analysis. Use disposable filter tips for adding reagents, DNA and control samples into tubes.

The total reaction volume is **25 µl**, the volume of the **cDNA** sample is **10 µl**.

1. Thaw the tubes with **PCR-mix-FL Influenza virus A H5, H7, H9**. Vortex the tubes with **PCR-mix-FL Influenza virus A H5, H7, H9, PCR-buffer-B** and **polymerase (TaqF)** and then centrifuge briefly.
2. Take the required number of tubes/strips for amplification of the cDNA obtained from test and control samples.
3. For N reactions, add to a new tube:

10*(N+1) µl of **PCR-mix-FL Influenza virus A H5, H7, H9**,

5*(N+1) µl of **PCR-buffer-B**

0.5*(N+1) µl of **polymerase (TaqF)** (see the scheme of reaction mixture preparation in Table 2).

4. Vortex the tube with prepared mixture, then centrifuge it briefly to sediment the drops.

Table 2

Scheme of reaction mixture preparation

Reagent volume per one reaction, μl	Reagent volume for specified number of reactions		
	10.0	5.0	0.5
Number of reactions*	PCR-mix-FL <i>Influenza virus A H5, H7, H9</i>	PCR-buffer-B	Polymerase (TaqF)
6	60	30	3.0
8	80	40	4.0
10	100	50	5.0
12	120	60	6.0
14	140	70	7.0
16	160	80	8.0
18	180	90	9.0
20	200	100	10.0
22	220	110	11.0
24	240	120	12.0
26	260	130	13.0
28	280	140	14.0
30	300	150	15.0
32	320	160	16.0

* Number of test samples including the controls of amplification, and one extra reaction (N+2+1).

5. Transfer **15 μl** of the prepared mixture to each tube.
6. Add **10 μl** of **cDNA samples** obtained at the RNA reverse transcription stage.
7. Carry out the control amplification reactions:

NCA	–	Add 10 μl of TE-buffer to the tube labeled NCA (Negative Control of Amplification).
C+	–	Add 10 μl of Positive Control <i>Influenza virus A H5, H7, H9</i> (C+<i>Influenza virus A H5, H7, H9</i>) to the tube labeled C+.

8.3.2 Amplification

1. Create a temperature profile on your instrument as follows:

Table 3

Amplification program

Step	Rotor-type instruments (e.g Rotor-Gene Q or equivalent)			Plate-type instruments (e.g CFX 96 Touch, CFX 96 Opus, QuantStudio 5 or equivalent.)		
	Temperature, °C	Time	Cycles	Temperature, °C	Time	Cycles
1	95	15 min	1	95	15 min	1
2	95	10 s	10	95	10 s	10
	54	20 s		54	25 s	
	72	10 s		72	25 s	
3	95	10 s	35	95	10 s	35
	54	20 s Fluorescence acquiring		54	25 s Fluorescence acquiring	
	72	10 s		72	25 s	

Fluorescent signal is detected in the channels for the **FAM, JOE** and **ROX** fluorophores.

2. Adjust the fluorescence channel sensitivity.
3. Insert tubes into the reaction module of the device.

NOTE: It is recommended to sediment drops from walls of tubes by short centrifugation (1–3 s) before placing them into the plate-type instrument.

4. Run the amplification program with fluorescence detection.
5. Analyze results after the amplification program is completed.

8.4 Instrument Settings

Test settings for rotor-type instruments

Channel	Calibrate/Gain Optimisation	Threshold	Dynamic tube	Slope Correct	More Settings/ Outlier Removal
FAM/Green	From 5FI to 10FI	0.1	On	-	10 %
JOE/Yellow	From 5FI to 10FI	0.1	On	-	10 %
ROX/Orange	From 5FI to 10FI	0.1	On	-	10 %

Test settings for plate-type instruments

Note: Set the heating/cooling rate **Ramp Rate** 2,5 °C/s by clicking the **Step Options** button for each step of cycling.

Channel	Threshold
FAM, JOE/HEX, ROX	Threshold is set at the level of 10-20 % of maximum fluorescence obtained for the Positive Controls in the last amplification cycle. Make sure that the fluorescence curve of the positive control has the typical exponential growth of fluorescence.

9 DATA ANALYSIS

Analysis of results is performed by the software of the real-time PCR instrument used by measuring fluorescence signal accumulation in three channels:

- The signal of amplification product of the haemagglutinin gene fragment cDNA of *Influenza virus A* subtype H5 is detected in the channel for the FAM fluorophore.
- The signal of amplification product of the haemagglutinin gene fragment cDNA of *Influenza virus A* subtype H7 is detected in the channel for the JOE fluorophore.
- The signal of amplification product of the haemagglutinin gene fragment cDNA of *Influenza virus A* subtype H9 is detected in the channel for the ROX fluorophore.

Results are interpreted by the crossing (or not-crossing) the fluorescence curve with the threshold line set at the specific level that corresponds to the presence (or absence) of a *Ct* value of the cDNA sample in the corresponding column of the results grid.

Principle of interpretation is the following:

- ***Influenza virus A* subtype H5 is identified** if the *Ct* value determined in the results grid for this sample in the channel for the FAM fluorophore is less than the specified boundary *Ct* value.
- ***Influenza virus A* subtype H7 is identified** if the *Ct* value determined in the results grid for this sample in the channel for the JOE fluorophore is less than the specified boundary *Ct* value.
- ***Influenza virus A* subtype H9 is identified** if the *Ct* value determined in the results grid for this sample in the channel for the ROX fluorophore is less than the specified boundary *Ct* value.
- The given *Influenza virus A* subtype is not identified (not detected) if the *Ct* values in the specified detection channel are absent.
- The result is **equivocal** if the *Ct* value determined in the respective channel is greater than the specified boundary *Ct* value. In this case, the PCR analysis of respective sample should be repeated. If the same result is obtained or the *Ct* value is determined less than threshold cycle, the sample is considered positive.

The result of the analysis is considered reliable only if the results obtained for Positive and Negative Controls of amplification are correct (see Table 4 and Table 5).

Table 4

Results for controls

Control	Stage for control	Ct value in the channel for fluorophore		
		FAM	JOE	ROX
NCA	PCR	Absent	Absent	Absent
C+	PCR	<boundary value	<boundary value	<boundary value
C-	RNA extraction	<boundary value	Absent	Absent

Table 5

Boundary Ct values

Sample	Rotor-type instruments			Plate-type instruments		
	Channel for fluorophore					
	FAM	JOE	ROX	FAM	JOE	ROX
NCA	Ct is absent			Ct is absent		
C+	<26	<26	<26	<26	<26	<26
Test samples	≤33	≤33	≤33	≤33	≤33	≤33

10 TROUBLESHOOTING

Results of analysis are not taken into account in the following cases:

1. If the Ct value determined for the Positive Control of Amplification (C+) in any channel is greater than the boundary Ct value or absent, the amplification should be repeated for all samples in which negative results was obtained in the respective channel.
2. If the Ct value is determined for the Negative Control of Amplification (NCA) in any channel, the amplification should be repeated for all samples in which positive result was obtained in the respective channel.

11 TRANSPORTATION

eSens Influenza virus A H5/H7/H9 QL PCR kit should be transported at 2–8 °C for no longer than 5 days.

12 STABILITY AND STORAGE

All components of the **eSens Influenza virus A H5/H7/H9 QL PCR kit** are to be stored at 2–8 °C when not in use (except for PCR-mix-FL *Influenza virus A H5, H7, H9*, PCR-buffer-B, and polymerase (TaqF)). All components of the **eSens Influenza virus A H5/H7/H9 QL PCR kit** are stable until the expiry date stated on the label. The shelf life of reagents before and after the first use is the same, unless otherwise stated.

NOTE: PCR-mix-FL *Influenza virus A H5, H7, H9*, PCR-buffer-B, and polymerase (TaqF) are to be stored at the temperature from minus 24 to minus 16 °C.

13 SPECIFICATIONS

13.1 Sensitivity

Clinical material	PCR kit	Sensitivity, GE/ml*
Nasal and oropharyngeal swabs, sputum (or tracheal aspirates), bronchoalveolar lavage, autopsy material in which the <i>Influenza virus A</i> RNA was detected	ES3303B	1 x 10 ³

* Genome equivalents (GE) of the pathogen agent per 1 ml of a sample transferred into the specified transport medium.

13.2 Specificity

The analytical specificity of **eSens Influenza virus A H5/H7/H9 QL PCR kit** is ensured by the selection of specific primers and probes as well as stringent reaction conditions. The primers and probes have been checked for possible homologies to all sequences published in gene banks by sequence comparison analysis.

The PCR kit detects the haemagglutinin genes fragments of the claimed *Influenza virus A* subtypes (H5, H7 and H9). The PCR kit specific activity is confirmed by analysis of strains of *Avian Influenza virus A/Anhui/1/2013* (H7N9), *A/Hong Kong/1073/99* (H9N2), *A/chicken/Moscow/2/07* (H5N1), and also by analysis of experimental samples of human material with addition of *Avian Influenza virus* strains.

The activity of PCR kit components was absent in respect of haemagglutinin genes fragments of *Influenza virus A* subtypes H1, H2, H3, H4, H13, H8, H6, H10, H11, H12, *Influenza virus B*, and also cDNA/DNA of strains and isolates of the main pathogens causing acute respiratory diseases as well as normal microflora of human nasal cavity and oropharynx and human DNA.

The clinical specificity of **eSens Influenza virus A H5/H7/H9 QL PCR kit** was confirmed in laboratory clinical trials.

13.3 Diagnostic characteristics

Results of PCR kit characteristics testing:

Samples description	Samples type	Number of samples	Results of using eSens Influenza virus A H5/H7/H9 QL PCR kit
Biological material containing <i>Influenza virus</i> A/H5 RNA*	Nasal and oropharyngeal swabs	100	Positive 100%
	Sputum	100	Positive 100%
	Autopsy material	100	Positive 100%
Biological material containing <i>Influenza virus</i> A/H7 RNA*	Nasal and oropharyngeal swabs	100	Positive 100%
	Sputum	100	Positive 100%
	Autopsy material	100	Positive 100%
Biological material containing <i>Influenza virus</i> A/H9 RNA*	Nasal and oropharyngeal swabs	100	Positive 100%
	Sputum	100	Positive 100%
	Autopsy material	100	Positive 100%
Biological material that does not contain <i>Influenza viruses</i> A/H5, A/H7, A/H9 RNA**	Nasal and oropharyngeal swabs	100	Negative 100%

* Model samples of biological material containing recombinant quality control samples was used as samples containing *Influenza viruses* A/H5, A/H7 and A/H9.

** Biological material samples from patients with suspected influenza containing *Influenza virus* A/H1N1pdm2009, *Parainfluenza viruses*, *Rhinoviruses* (that was proved by testing for example with **eSens ARVI screen QL PCR kit**) was used as samples that do not contain *Influenza viruses* A/H5, A/H7 and A/H9.

In accordance with the submitted data the **diagnostic sensitivity** of the **eSens Influenza virus A H5/H7/H9 QL PCR kit** is 98-100 % with a confidence coefficient of 90 % for all type of the biological material.

The **diagnostic specificity** of the **eSens Influenza virus A H5/H7/H9 QL PCR kit** is 98-100 % with a confidence coefficient of 90 %.

14 QUALITY CONTROL

The production process, including batch release, is carried out in accordance with an established quality management system certified according to ISO 13485.

15 KEY TO SYMBOLS USED

	Catalogue number		Caution
	Batch code		Contains sufficient for <n> tests
	<i>In vitro</i> diagnostic medical device		Use-by Date
	Version		Consult instructions for use
	Temperature limit		Keep away from sunlight
	Manufacturer	NCA	Negative control of amplification
	Date of manufacture	C+	Positive control of amplification
	Authorized representative in the European Community		

List of Changes Made in the Instruction Manual

VER	Location of changes	Essence of changes
01_04/2022		

Ecoli Dx, s.r.o. , Purkyňova 74/2



110 00 Praha 1, Česká republika
Tel: +420 325 209 912

Mobil: +420 739 802 523