

MAGNO-sorb Nucleic Acid Extraction Kit

Instruction Manual

KEY TO SYMBOLS USED

| | | | |
|---|---|---|-----------------------------------|
|  | Catalogue number |  | Contains sufficient for <n> tests |
|  | Batch code |  | Use-by Date |
|  | In vitro diagnostic medical device |  | Consult instructions for use |
|  | Version |  | Negative Control of extraction |
|  | Temperature limit |  | Positive Control of extraction |
|  | Manufacturer |  | Internal Control |
|  | Date of manufacture |  | GHS02: Flame |
|  | Authorized representative in the European Community |  | GHS05: Corrosion |
|  | Caution |  | GHS07: Exclamation mark |

1. INTENDED USE

MAGNO-sorb nucleic acid extraction kit is intended for extraction of RNA/DNA from the human biological material for subsequent testing by the nucleic acid amplification techniques (NAT):

- blood plasma,
- whole blood,
- cerebrospinal fluid,
- discharge of the conjunctiva,
- saliva,
- nasopharyngeal and oropharyngeal swabs,
- sputum / endotracheal aspirate,
- bronchoalveolar lavage fluid / bronchial washings,
- pleural / ascitic fluid,
- pus / liquid or dense necrotic content,
- tissue (biopsy, surgical, autopsy) native material,
- paraffin-embedded tissue (biopsy, surgical, autopsy) material,
- urogenital mucous discharge (swab, scrape) (vagina, cervical canal (ectocervix and endocervix), urethra),
- urine,
- feces,
- anal canal/rectal mucosal discharge (swab, scrape).

The reagent kit can be used for extraction of total RNA/DNA from ticks, food products, water (drinking water, open sources, sewage), germ cultures, washings from environmental objects for sanitary and anti-epidemic (preventive) measures.

Indications and contra-indications for use of the reagent kit

DNA/RNA extraction is used in preanalytical stage of in vitro diagnostics by NAT.

2. PRINCIPLE OF NUCLEIC ACID EXTRACTION

A test sample is treated by the lysis solution in the presence of the magnetic silica particles (magnetic sorbent). As a result the cell membranes, viral envelopes and other biopolymer complexes are destructed and the nucleic acids are released. The dissolved nucleic acids bind to the sorbent particles while other components of the lysed biological material stay in the solution and are removed by sorbent precipitation at magnetic rack and subsequent washings. The nucleic acids are transferred from the sorbent surface to the solution after adding the buffer for elution to the magnetic sorbent. Then the solution is separated from the sorbent by magnetic force.

The obtained nucleic acid sample is highly purified and free from inhibitors of amplification, which provides high analytical sensitivity of NAT assay.

3. CONTENT

MAGNO-sorb nucleic acid extraction kit is produced in 4 forms:

- variant 100-200,  K3-1061-100-CE,
- variant 100-1000,  K3-1062-100-CE,
- variant 100-100M,  K3-1063-100-CE,
- variant 100-200M,  K3-1064-100-CE.

Variant 100-200 or 100-1000 includes:

| Reagent | Description | variant 100-200 | | variant 100-1000 | |
|---------------------------|---|-----------------|----------|------------------|----------|
| | | Volume, ml | Quantity | Volume, ml | Quantity |
| Lysis Solution MAGNO-sorb | clear liquid from colorless to yellow or pink colour ¹ | 107 | 1 vial | 90 | 3 vials |
| Component A | colorless clear liquid | 1.0 | 1 tube | 1.0 | 2 tubes |
| Washing Solution 5 | clear liquid from colorless to yellow or pink colour ¹ | 75 | 2 vials | 75 | 2 vials |
| Washing Solution 6 | colorless clear liquid | 87 | 1 vial | 87 | 1 vial |
| Washing Solution 7 | colorless clear liquid ² | 35 | 1 vial | 35 | 1 vial |
| Magnetized silica | suspension of magnetic particles | 1.0 | 2 tubes | 1.5 | 2 tubes |
| Buffer for elution A | colorless clear liquid | 1.25 | 8 tubes | 1.25 | 8 tubes |

Variant 100-200 is intended for DNA/RNA extraction from 100 samples (including controls). The volume of test material is 200 µl.

Variant 100-1000 is intended for DNA/RNA extraction from 100 samples (including controls). The volume of test material is 1,000 µl.

Variant 100-100M or 100-200M includes:

| Reagent | Description | variant 100-100M | | variant 100-200M | |
|---------------------------|---|------------------|----------|------------------|----------|
| | | Volume, ml | Quantity | Volume, ml | Quantity |
| Lysis Solution MAGNO-sorb | clear liquid from colorless to yellow or pink colour ¹ | 40 | 1 vial | 70 | 1 vial |
| Component A | colorless clear liquid | 0.5 | 1 tube | 1.0 | 1 tube |
| Washing Solution 5 | clear liquid from colorless to yellow or pink colour ¹ | 60 | 1 vial | 60 | 1 vial |
| Washing Solution 4 | colorless clear liquid | 60 | 1 vial | 60 | 1 vial |
| Magnetized silica | suspension of magnetic particles | 1.0 | 1 tube | 1.0 | 2 tubes |
| Buffer for elution A | colorless clear liquid | 1.25 | 8 tubes | 1.25 | 8 tubes |

Variant 100-100M is intended for DNA/RNA extraction from 100 samples (including controls). The volume of test material is 100 µl.

Variant 100-200M is intended for DNA/RNA extraction from 100 samples (including controls). The volume of test material is 200 µl.

4. ADDITIONAL REQUIREMENTS

- 1.5-ml disposable polypropylene screwed or tightly closed tubes.
- 5-ml disposable polypropylene or polystyrene round-bottomed tubes, 12 mm diameter - when using variant 100-1000.
- 5-ml disposable polypropylene or polystyrene tubes, 12 mm diameter - when using variant 100-1000.
- Sterile DNase- and RNase-free pipette tips without filter (up to 200, 1,000, and 5,000 µl).
- Sterile DNase- and RNase-free pipette filter tips (from 200 to 1,000 µl, and 5,000 µl).
- Tube racks.
- Magnetic racks for 1.5-ml tubes.
- Magnetic racks for 5-ml tubes (12 mm diameter) - when using variant 100-1000.
- Pipettes (adjustable).
- Desktop microcentrifuge.
- PCR box.
- Vortex mixer.
- Thermostat for 1.5-ml tubes with capable of incubating at 25-100 °C.
- Thermostat for 5-ml tubes (12 mm diameter) with capable of incubating at 25-100 °C - when using variant 100-1000.
- Vacuum aspirator with flask for removing supernatant.
- Refrigerator at the temperature from 2 to 8 °C.
- Disposable powder-free gloves and a laboratory coat.
- Reservoir for used tips.

In case of using automated stations for nucleic acid extraction:

- Open automated DNA extraction systems with magnetic stirring (for example, Auto-Pure 96, Auto-Pure 24 (Hangzhou Allsheng Instruments Co., Ltd., China), KingFisher Flex (Thermo Fisher Scientific Oy), Nexor 96 (Yantai Addcare Bio-Tech Co., Ltd., China)).
- Open automated DNA extraction pipetting systems (for example, Xiril (Neon 100 series) (Tecan Schweiz AG, Switzerland), MicroLab STARlet (Hamilton Bonaduz AG, Switzerland)).
- Set of consumables for used automated system according to the manufacturer's recommendations including consumables for stations with magnetic stirring:
 - 96-tip comb - when using variant 100-200, variant 100-100M, variant 100-200M or 24-tip comb - when using variant 100-1000;
 - 96-deepwell plate - when using variant 100-200, variant 100-100M, variant 100-200M or 24-deepwell plate - when using variant 100-1000;
 - 96-elution plate - when using variant 100-200, variant 100-100M, variant 100-200M; - sterile self-adhesive film for plates.
- 1.5-5 ml disposable polypropylene or polystyrene tubes with caps - to prepare the required volume of a mixture of Magnetized silica.
- Sterile pipette filter tips (up to 100 µl).
- Sterile pipette filter tips (up to 1,000 and 5,000 µl).
- Single-channel pipettes (adjustable).
- 8-channel pipettes.
- Reservoirs (baths) for filling 8-channel dispensers.

¹ If Lysis Solution MAGNO-sorb and Washing Solution 5 are stored below 20 °C, crystalline precipitate may form.

² An empty vial is provided. Before the work with nucleic acid extraction kit, add 35 ml of acetone (mass fraction of acetone should be no less than 99.7 %) into the vial. Work in a fume hood.

WARNING! In order to prevent acetone evaporation the vial cap should be tightly closed.

5. GENERAL PRECAUTIONS

The user should always pay attention to the following:

- Use sterile DNase- and RNase-free pipette filter tips and use new tip for every procedure.
- Use disposable protective gloves and laboratory cloths, and protect eyes while samples and reagents handling. Thoroughly wash hands afterwards.
- Do not eat, drink, smoke, apply cosmetics, or handle contact lenses in laboratory work areas.
- Do not use the kit if the internal packaging was damaged or its appearance was changed.
- Do not use the kit if the transportation and storage conditions according to the Instruction Manual were not observed.
- Do not use a kit after its expiration date.
- Dispose of all specimens and unused reagents in accordance with local regulations.
- Samples should be considered potentially infectious and handled in a biological cabinet in compliance with appropriate biosafety practices.
- Clean and disinfect all sample or reagent spills using a disinfectant, such as 0.5% sodium hypochlorite, or another suitable disinfectant.
- Avoid inhalation of vapors, samples and reagents contact with the skin, eyes, and mucous membranes. Harmful if swallowed. If these solutions come into contact, rinse the injured area immediately with water and seek medical advice if necessary.
- While observing the conditions of transportation, operation and storage, there are no risks of explosion and ignition.
- Safety Data Sheets (SDS) are available on request.
- The kit is intended for single use for extraction of specified number of samples (see the section "Content").
- The kit is ready for use in accordance with the Instruction Manual. Use the PCR kit strictly for intended purpose.
- Use of this product should be limited to personnel trained in the techniques of DNA/RNA extraction.
- Workflow in the laboratory must be one-directional, beginning in the Extraction Area and moving to the Amplification and Detection Area. Do not return samples, equipment and reagents to the area where the previous step was performed.



Some components of this kit contain sodium azide as a preservative. Do not use metal tubing for reagent transfer.

| | |
|---|---|
| <p>Lysis Solution MAGNO-sorb</p> <p>Danger</p> | <p>Contains substance: isopropanol, guanidine thiocyanate, Triton X-100, 1-Thioglycerol</p> <p>H226: Flammable liquid and vapour H302: Harmful if swallowed. H314: Causes severe skin burns and eye damage. H318: Causes serious eye damage. H332: Harmful if inhaled. H336: May cause drowsiness or dizziness. H412: Harmful to aquatic life with long lasting effects.</p> <p>EUH032: Contact with acids liberates very toxic gas.</p> <p>P210: Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. P260: Do not breathe vapours. P273: Avoid release to the environment. P305+P351+P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses if present and easy to do – continue rinsing. P501: Dispose of contents in accordance with national regulation.</p> |
| <p>Washing Solution 5</p> <p>Danger</p> | <p>Contains substance: isopropanol, guanidine thiocyanate.</p> <p>H226: Flammable liquid and vapour H302: Harmful if swallowed. H314: Causes severe skin burns and eye damage. H318: Causes serious eye damage. H332: Harmful if inhaled. H336: May cause drowsiness or dizziness H412: Harmful to aquatic life with long lasting effects.</p> <p>EUH032: Contact with acids liberates very toxic gas.</p> <p>P210: Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. P260: Do not breathe vapours. P280: Wear protective gloves/protective clothing/eye protection/face protection. P305+P351+P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses if present and easy to do – continue rinsing. P303+P361+P353: IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water. P501: Dispose of contents in accordance with national regulation.</p> |
| <p>Washing Solution 4</p> <p>Danger</p> | <p>Isopropanol EC No 200-661-7 CAS No 67-63-0</p> <p>H225: Highly flammable liquid and vapour. H319: Causes serious eye irritation. H336: May cause drowsiness or dizziness.</p> <p>P210: Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. P261: Avoid breathing vapours. P280: Wear protective gloves/protective clothing/eye protection/face protection. P305+P351+P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses if present and easy to do – continue rinsing. P403+P233: Store in a well ventilated place. Keep container tightly closed. P501: Dispose of contents in accordance with national regulation.</p> |
| <p>Washing Solution 6</p> <p>Warning</p> | <p>Contains substance: Isopropanol</p> <p>H226: Flammable liquid and vapour. H319: Causes serious eye irritation. H336: May cause drowsiness or dizziness.</p> <p>P210: Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. P261: Avoid breathing vapours. P264: Wash your hand thoroughly after handling. P305+P351+P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses if present and easy to do – continue rinsing. P403+P233: Store in a well ventilated place. Keep container tightly closed. P501: Dispose of contents in accordance with national regulation.</p> |
| <p>Washing Solution 7</p> | <p>Due to providing the empty glass vial, the labeling of Washing Solution 7 as a hazard reagent is absent. Then working with acetone, see hazard and precautionary statements provided by a supplier.</p> |

6. SAMPLING AND HANDLING

See the information about the sampling, conditions of transportation and storage

NOTE: of the test material, as well as the necessity and procedure of its pretreatment before DNA/RNA extraction in the Instruction manual for the PCR kit.

MAGNO-sorb nucleic acid extraction kit is recommended for DNA and RNA extraction from:

- blood plasma,
- whole blood,
- cerebrospinal fluid,
- discharge of the conjunctiva,
- saliva,
- nasopharyngeal and oropharyngeal swabs,
- sputum / endotracheal aspirate,
- bronchoalveolar lavage fluid / bronchial washings,
- pleural / ascitic fluid,
- pus / liquid or dense necrotic content,
- tissue (biopsy, surgical, autopsy) native material,
- paraffin-embedded tissue (biopsy, surgical, autopsy) material,
- urogenital mucous discharge (swab, scrape) (vagina, cervical canal (ectocervix and endocervix), urethra),
- urine,
- feces,
- anal canal/rectal mucosal discharge (swab, scrape).

The reagent kit can be used for extraction of total RNA/DNA from ticks, food products, water (drinking water, open sources, sewage), germ cultures, washings from environmental objects for sanitary and anti-epidemic (preventive) measures.

Interfering substances and limitations of using test material samples

The information about limitations of using test material samples is specified in the Instruction Manual of the PCR kit.

Potential interfering substances

Endogenous and/or exogenous substances that may be present in the biological material used for the study were selected to assess potential interference (see Table 1).

Model samples of various biological material without adding and with the addition of potentially interfering substances were tested. The maximum concentration of potentially interfering substances in model samples and interference presence are listed in Table 1.

| Test material | Type of potential interferent | Potential interferent | Tested concentration in a sample | Interference presence |
|---|-------------------------------|--|----------------------------------|-----------------------|
| Saliva, nasopharyngeal and oropharyngeal swabs, urogenital mucous discharge (swab, scrape) (vagina, cervical canal (ectocervix and endocervix), urethra) | Exogenous substances | Chlorhexidine bigluconate aqueous solution | 2.5% | Not detected |
| Saliva, nasopharyngeal and oropharyngeal swabs, sputum/ endotracheal aspirate, bronchoalveolar lavage fluid/bronchial washings, anal canal/rectal mucosal discharge (swab, scrape), urogenital mucous discharge (swab, scrape) (vagina, cervical canal (ectocervix and endocervix), urethra) | Endogenous substances | Mucin | 9 mg/ml | Not detected |
| Blood plasma, saliva, nasopharyngeal and oropharyngeal swabs, bronchoalveolar lavage fluid/bronchial washings, tissue (autopsy) native material, anal canal/rectal mucosal discharge (swab, scrape), urogenital mucous discharge (swab, scrape) (vagina, cervical canal (ectocervix and endocervix), urethra) | Endogenous substances | Hemoglobin | 0.21 g/l | Not detected |
| Sputum/ endotracheal aspirate, bronchoalveolar lavage fluid / bronchial washings, urine, pleural / ascitic fluid, pus / liquid or dense necrotic contents, feces, cerebrospinal fluid | Exogenous substances | "Rifampicin" | 9-17.5 µg/ml | Not detected |
| | | "Isoniazid" | 3-7 µg/ml | Not detected |
| | | "Ethambutol" | 2-5 µg/ml | Not detected |
| | | "Pyrazinamide" | 39 µg/ml | Not detected |
| | | "Ofloxacin" | 5.2 µg/ml | Not detected |
| | | "Kanamycin" | 22 µg/ml | Not detected |
| | | Para-aminosalicylic acid (PAS) | 75 µg/ml | Not detected |
| Endogenous substances | Blood | sample : blood 1:1 | Not detected | |
| Paraffin-embedded tissue (biopsy, surgical, autopsy) material | Exogenous substances | Formalin sour | — | <u>Detected</u> |
| | | Buffered formalin | — | Not detected |
| Whole blood | Endogenous substances | Hemoglobin | 250 g/l | Not detected |
| | | Total bilirubin | 210 µmol/l | Not detected |
| | | Total cholesterol | 78 mmol/l | Not detected |
| | | Triglycerides | 37.0 mmol/l | Not detected |
| Discharge of the conjunctiva | Exogenous substances | Lithium heparin | 12-30 IU/ml | <u>Detected</u> |
| | | Potassium EDTA | 2.0 µg/ml | Not detected |
| Discharge of the conjunctiva | Exogenous substances | "Dicain®" | 0.003% | Not detected |

| Test material | Type of potential interferent | Potential interferent | Tested concentration in a sample | Interference presence |
|--|-------------------------------|--|---|-----------------------|
| Urogenital mucous discharge (swab, scrape) (vagina, cervical canal (ectocervix and endocervix), urethra) | Endogenous substances | Lactoferrin | 5 µg/ml | Not detected |
| | | Glycogen | 120 mg/ml | Not detected |
| | Exogenous substances | "Clotrimazole" | 16 % (the drug volume to the volume of the test sample) | Not detected |
| | | "Metronidazole" | | Not detected |
| | | "Miramistin [®] " | | Not detected |
| | | "Polygynax" | | Not detected |
| | | "Macmiror [®] Complex" | | Not detected |
| | | "Contex Silk", silicone gel lubricant for intimate use | | Not detected |

7. WORKING CONDITIONS

MAGNO-sorb nucleic acid extraction kit should be used at the temperature from 20 to 28 °C and relative humidity from 15 to 75 %.

8. PROTOCOL

8.1 Variant 100-200

The volume of test sample for the extraction is 200 µl.

Test material:

- blood plasma,
- whole blood,
- cerebrospinal fluid,
- nasopharyngeal and oropharyngeal swabs,
- sputum / endotracheal aspirate,
- bronchoalveolar lavage fluid / bronchial washings,
- anal canal/rectal mucosal discharge (swab, scrape),
- urine,
- blood-sucking arthropods (ticks),
- water (drinking water, open sources, sewage).

NOTE: See the information about pretreatment of biological material for DNA/RNA extraction in the Instruction manual for the used PCR kit.

NOTE: To work with RNA, it is necessary to use only disposable sterile plastic consumables that have a special marking RNase-free, DNase-free.

NOTE: **Lysis Solution MAGNO-sorb** has unpleasant odor. Work in the PCR box.

Manual procedure when using magnetic racks

1. If there is a precipitate in **Lysis Solution MAGNO-sorb** and **Washing Solution 5** warm them up at 60 °C until crystals disappear. Resuspend **Magnetized silica** on vortex.
2. Prepare the required number of disposable 1.5-ml tubes (including a tube for the Negative Control of extraction (C) and a tube for the Positive Control of Extraction (PCE) if they are provided for analysis). Mark the tubes.
3. Mix in a separate tube **Internal Control (IC)** (if it is provided for analysis), **Component A**, and **Magnetized silica** in the following proportion calculated per one sample: **10 µl of Internal Control (IC)**, **10 µl of Component A**, and **20 µl of Magnetized Silica**. Do not forget to add extra volumes for one more reaction. For example:

| Number of samples for DNA/RNA extraction | Internal Control (IC), µl | Component A, µl | Magnetized silica, µl |
|--|---------------------------|-----------------|-----------------------|
| 6 | 70 | 70 | 140 |
| 12 | 130 | 130 | 260 |
| 18 | 190 | 190 | 380 |
| 24 | 250 | 250 | 500 |

NOTE: It is allowed to change the volume of **Internal Control (IC)** according to the Instruction manual for the used PCR kit.

4. Add **40 µl** of the resuspended mixture of **Internal Control (IC)**, **Component A**, and **Magnetized silica** into each tube (recalculate the volume of the mixture if the volume of Internal Control (IC) has been changed);

5. Add **900 µl** of **Lysis Solution MAGNO-sorb** into the tubes.

6. Add **200 µl** of **test sample** into each tube, using a separate filter tip for each sample. Mix by vortexing.

NOTE: It is allowed to change the test sample volume according to the Instruction manual for the used PCR kit.

7. Add **200 µl** of the **Negative Control (C-)** reagent into the tube for the Negative Control of extraction (C-) and **200 µl** of **required Positive Control sample** into the tube for the Positive Control of extraction (PCE) (if they are provided for PCR). Tightly close the tubes and thoroughly mix by vortexing.

NOTE: It is allowed to change the volume and dilution of **Negative Control (C-)** and **Positive Control** according to the Instruction manual for the used PCR kit.

8. Incubate the tubes at **60 °C** for **10 min**.

9. Centrifuge the tubes shortly, transfer the tubes to a magnetic rack, and incubate for **2 min**.

NOTE: If it is inconvenient/impossible to open the tubes in a magnetic rack without sorbent resuspension, it is necessary to open the tubes before placing them in a magnetic rack.

10. Leaving the tubes in the magnetic rack, carefully remove the supernatant inserting a tip near the internal tube wall and using vacuum aspirator. Take a new tip for each sample.

11. Add **700 µl** of **Washing Solution 5** to the tubes. Tightly close the tubes.

NOTE: If it is inconvenient/impossible to close tightly the tubes in a magnetic rack, it is necessary to place the tubes in a regular tube rack.

12. Resuspend **Magnetized silica** by mixing on vortex. Then sediment the drops by short centrifugation.

13. Place the tubes in a magnetic rack for **2 min**.

14. Open the tubes, then carefully remove the supernatant as described above.

15. Repeat washing procedure with **700 µl** of **Washing Solution 5** (steps 11-14).

16. Carry out washing procedure with **700 µl** of **Washing Solution 6** as described above.

17. Add **200 µl** of **Washing Solution 7** to each tube, mix, and vortex shortly to sediment drops.

18. Transfer the tubes into the magnetic rack for **1 min**, open the tubes and remove the supernatant.

19. Dry the magnetized silica. To do this, leave the tubes with open caps in the magnetic rack for **10-20 min**.

20. Add **100 µl** of **Buffer for elution A** to each tube, close the tubes and mix on vortex.

NOTE: It is allowed to change the elution volume according to the Instruction manual for the used PCR kit.

21. Incubate the tubes at **60 °C** for **5 min**. Vortex the tubes in 2 min.

22. Sediment the drops by short centrifugation and transfer the tubes to the magnetic rack. Incubate for **2 min**. Supernatant contains purified DNA and RNA.

NOTE: Do not take the tubes away from a magnetic rack when removing DNA/RNA.

Procedure for automated stations for nucleic acid extraction

The detailed information work with automated stations see in the *Instruction Manual* given by the manufacturer of the instrument.

Extraction using stations with magnetic stirring (for example, Auto-Pure 96 (Hangzhou Allsheng Instruments Co., Ltd., China), KingFisher Flex (Thermo Fisher Scientific Oy, Finland)).

Only the protocols located on the amplisens.ru website should be used for extraction.

NOTE: The extraction protocol for the automated station contains specific parameters that are not available for manual extraction.

1. If there is a precipitate in **Lysis Solution MAGNO-sorb** and **Washing Solution 5** warm them up at 60 °C until crystals disappear. Resuspend **Magnetized silica** on vortex.
2. Prepare the required number of 96-well plates (see point 3) compatible with the instrument; mark them according to the protocol.

3. Add in the plates according to the marking and quantity of samples for extraction:
 - **760 µl** of **Lysis Solution MAGNO-sorb** in a deepwell plate,
 - **700 µl** of **Washing Solution 5** per two deepwell plates,
 - **700 µl** of **Washing Solution 6** in a deepwell plate,
 - **200 µl** of **Washing Solution 7** in a deepwell plate,
 - **100 µl** of **Buffer for elution A** in the plate for elution.

NOTE: It is allowed to change the elution volume according to the Instruction manual for the used PCR kit

4. Mix in a separate tube in the following proportion calculated per one sample:
 - **10 µl** of **Internal Control (IC)** (if it is provided for analysis),
 - **10 µl** of **Component A**,
 - **20 µl** of **Magnetized silica**.

Prepare mixture for total number of test and control samples plus extra samples (for example, in case of extraction from 96 samples it is recommended to prepare mixture for 100 extractions, i.e. plus 4 extra samples).

NOTE: It is allowed to change the volume of **Internal Control (IC)** according to the Instruction manual for the used PCR kit.

5. Add **40 µl** of thoroughly resuspended mixture to each well with **Lysis Solution MAGNO-sorb** (recalculate the volume of the mixture if the volume of Internal Control (IC) has been changed).

6. Add **200 µl** of **test and control** (if they are provided for analysis) **samples** in each well with **Lysis Solution MAGNO-sorb** and mixture according to the marking, using a separate filter tip for each sample.

NOTE: It is allowed to change test and control sample volumes according to the Instruction manual for the used PCR kit

7. Choose the protocol of extraction, place the prepared plates and tip-plate according to the protocol and placement indicators.

8. Run the nucleic acid extraction protocol.

9. If necessary, seal the elution plate with the purified nucleic acid using a self-adhesive film at the end of the extraction.

Extraction using pipetting stations (for example, Xiril (Neon 100 series) (Tecan Schweiz AG, Switzerland), MicroLab STARlet (Hamilton Bonaduz AG, Switzerland)).

1. Arrange the required number of the consumables on the desktop of the station: tips, test tubes/plates for elution.

2. If there is a precipitate in **Lysis Solution MAGNO-sorb** and **Washing Solution 5** warm them up at 60 °C until crystals disappear. Resuspend **Magnetized silica** on vortex.

3. Place the tubes with **Internal Control (IC)** (if it is provided for analysis), **Component A**, **Buffer for elution A** and **Magnetized silica** on the desktop of the station.

4. Pour **Lysis Solution MAGNO-sorb**, **Washing Solution 5**, **Washing Solution 6** and **Washing Solution 7** into the appropriate containers placed on the desktop of the station.

5. Place the tubes with test samples on the desktop of the station. Choose the rack according to the type of used tubes.

6. Add no less than **150 µl** of **Negative Control (C-)** into the tube for the Negative Control of extraction (C-) and no less than **150 µl** of **required Positive Control sample** into the tube for the Positive Control of extraction (PCE) (if they are provided for PCR).

NOTE: It is allowed to change the volume and dilution of **Negative Control (C-)** and **Positive Control sample** according to the Instruction manual for the used PCR kit. The volume of control samples should be increased taking into account some extra volume.

7. Choose and run the protocol of extraction.

8. If the elution was carried out into a plate, seal the elution plate with the purified nucleic acid using a self-adhesive film at the end of the extraction if necessary.

Storage of purified nucleic acids

Purified DNA samples can be stored:

- at the temperature from 2 to 8 °C for 1 week;

- at the temperature from minus 24 to minus 16 °C for 1 year.

Purified RNA samples can be stored:

- at the temperature from 2 to 8 °C for 4 hours;

- at the temperature from minus 24 to minus 16 °C for 1 month;

- at the temperature not more than minus 68 °C for 1 year.

In case of manual extraction transfer the supernatant for purified nucleic acids storage into a new tube without disturbing the sediment.

In case of automatic extraction transfer eluate into the tubes for purified nucleic acids storage at subzero temperature.

8.2 Variant 100-1000

The volume of test sample for the extraction is 1,000 µl.

Test material:

- blood plasma,
- cerebrospinal fluid,
- urine,
- water (drinking water, open sources, sewage).

NOTE: See the information about pretreatment of biological material for DNA/RNA extraction in the Instruction manual for the used PCR kit.

NOTE: To work with RNA, it is necessary to use only disposable sterile plastic consumables that have a special marking RNase-free, DNase-free.

NOTE: **Lysis Solution MAGNO-sorb** has unpleasant odor. Work in the PCR box.

Manual procedure when using magnetic racks

- If there is a precipitate in **Lysis Solution MAGNO-sorb** and **Washing Solution 5** warm them up at 60 °C until crystals disappear. Resuspend **Magnetized silica** on vortex.
- Prepare the required number of disposable 5-ml tubes (including a tube for the Negative Control of extraction (C) and a tube for the Positive Control of Extraction (PCE) if provided with the amplification kit) and tube caps. Mark the tubes.
- Mix in a separate tube **Internal Control (IC)** (if it is provided for analysis), **Component A**, and **Magnetized silica** in the following proportion calculated per one sample: **10 µl of Internal Control IC**, **20 µl of Component A**, and **30 µl of Magnetized Silica**. Do not forget to add extra volumes for one more reaction. For example:

| Number of samples for DNA/RNA extraction | Internal Control (IC), µl | Component A, µl | Magnetized silica, µl |
|--|---------------------------|-----------------|-----------------------|
| 6 | 70 | 140 | 210 |
| 12 | 130 | 260 | 390 |
| 18 | 190 | 380 | 570 |
| 24 | 250 | 500 | 750 |

NOTE: It is allowed to change the volume and dilution of **Internal Control (IC)** according to the Instruction manual for the used PCR kit.

- Add **60 µl** of the thoroughly resuspended mixture of **Internal Control (IC)**, **Component A**, and **Magnetized silica** into each 5-ml tube (recalculate the volume of the mixture if the volume of Internal Control (IC) has been changed);
- Add **2.6 ml of Lysis Solution MAGNO-sorb** into the tubes.
- Add **1 ml of test sample** to each tube using a separate filter tip for each sample. Mix by pipetting. Cap the tubes.

NOTE: It is allowed to change the test sample volume according to the Instruction manual for the used PCR kit.

- Add **1 ml of the Negative Control (C-)** into the tube for the Negative Control of extraction (C-) and **1 ml of required Positive Control sample** into the tube for the Positive Control of extraction (PCE) (if they are provided for PCR). Mix by pipetting. Cap the tube.

NOTE: It is allowed to change the volume and dilution of **Negative Control (C-)** and **Positive Control** according to the Instruction manual for the used PCR kit.

- Incubate the tubes at **60 °C for 10 min**.
- Open the tubes; transfer them to a magnetic rack and incubate for **6 min**.
- Leaving the tubes in the magnetic rack, carefully remove the supernatant inserting the tip near the internal tube wall and using vacuum aspirator. Take a new 1,000-µl tip for each sample. Transfer the tubes to a regular tube rack.
- Add **700 µl of Washing Solution 5** to the tubes. Cap the tubes.
- Take the required number of disposable 1.5-ml tubes (including the tubes for the Positive and Negative Controls of extraction). Mark the tubes.
- Vortex the tubes and then pipette to remove magnetic beads from tube walls. Transfer the entire content of the tubes to the prepared 1.5-ml tubes.
- Place the 1.5-ml tubes in the magnetic rack and incubate for 2 min.

NOTE: If it is inconvenient/impossible to open the tubes in a magnetic rack without sorbent resuspension, it is necessary to open the tubes before placing them in a magnetic rack.

- Carefully remove the supernatant as described above.
- Add **700 µl of Washing Solution 5**, close the tubes.

NOTE: If it is inconvenient/impossible to close tightly the tubes in a magnetic rack, it is necessary to place the tubes in a regular tube rack.

- Resuspend **Magnetized silica** by mixing on vortex. Then sediment the drops by short centrifugation and then repeat steps 14-15.
- Carry out washing procedure with **700 µl of Washing Solution 6** as described above.
- Add **200 µl of Washing Solution 7**, mix, and vortex shortly to remove drops.
- Place the tubes to the magnetic rack for 1 min and then carefully remove the supernatant.
- Dry the magnetized silica. To do this, open the tubes and incubate them in the magnetic rack for 10-20 min.
- Add **100 µl of Buffer for elution A** to each tube, close the tubes and mix on vortex.

NOTE: It is allowed to change the elution volume according to the Instruction manual for the used PCR kit.

- Incubate the tubes at **60 °C for 5 min**. Vortex in 2 min.
 - Sediment the drops by short centrifugation and transfer the tubes to the magnetic rack. Incubate for 2 min.
- Supernatant contains purified DNA and RNA.

NOTE: Do not take the tubes away from a magnetic rack when removing DNA/RNA.

Procedure for automated stations for nucleic acid extraction

The detailed information work with automated stations see in the Instruction Manual given by the manufacturer of the instrument.

Extraction using stations with magnetic stirring (for example, KingFisher Flex (Thermo Fisher Scientific Oy, Finland), Auto-Pure 24 (Hangzhou Allsheng Instruments Co., Ltd., China), Nexor 96 (Yantai Addcare Bio-Tech Co., Ltd., China).

Only the protocols located on the amplicons.ru website should be used for extraction.

NOTE: The extraction protocol for the automated station contains specific parameters that are not available for manual extraction.

- If there is a precipitate in **Lysis Solution MAGNO-sorb** and **Washing Solution 5** warm them up at 60 °C until crystals disappear. Resuspend **Magnetized silica** on vortex.
- Prepare the required number of 24-well plates (see point 3) compatible with the instrument; mark them according to the protocol.
- Add in the plates according to the marking and quantity of samples for extraction:
 - 2.6 ml of Lysis Solution MAGNO-sorb** in a deepwell plate,
 - 700 µl of Washing Solution 5** per two deepwell plate,
 - 700 µl of Washing Solution 6** in a deepwell plate,
 - 200 µl of Washing Solution 7** in a deepwell plate,
 - 100 µl of Buffer for elution A** in the plate for elution.

NOTE: It is allowed to change the elution volume according to the Instruction manual for the used PCR kit.

- Mix in a separate tube in the following proportion calculated per one sample:
 - 10 µl of Internal Control (IC)** (if it is provided for analysis),
 - 20 µl of Component A**,
 - 30 µl of Magnetized silica**.

Prepare mixture for total number of test and control samples plus extra samples (for example, in case of extraction from 24 samples it is recommended to prepare mixture for 26 extractions, i.e. plus 2 extra samples).

NOTE: It is allowed to change the volume and dilution of **Internal Control (IC)** according to the Instruction manual for the used PCR kit.

- Add **60 µl** of thoroughly resuspended mixture to each well with **Lysis Solution MAGNO-sorb** (recalculate the volume of mixture if the volume of Internal Control (IC) has been changed).
- Add **1 ml of test and control** (if they are provided for analysis) **samples** in each well with **Lysis Solution MAGNO-sorb** and mixture according to the marking using a separate filter tip for each sample.

NOTE: It is allowed to change test and control sample volumes according to the Instruction manual for the used PCR kit.

- Choose the protocol of extraction, place the prepared plates and tip-plate according to the protocol and placement indicators.
- Run the nucleic acid extraction protocol.
- If necessary, seal the elution plate with the purified nucleic acid using a self-adhesive film at the end of the extraction.

Extraction using pipetting stations (for example, Xiril (Neon 100 series) (Tecan Schweiz AG, Switzerland), MicroLab STARlet (Hamilton Bonaduz AG, Switzerland)).

- Arrange the required number of the consumables on the desktop of the station: tips, test tubes/elution plates.
- If there is a precipitate in **Lysis Solution MAGNO-sorb** and **Washing Solution 5** warm them up at 60 °C until crystals disappear. Resuspend **Magnetized silica** on vortex.
- Place the tubes with **Internal Control (IC)** (if it is provided for analysis), **Component A**, **Buffer for elution A** and **Magnetized silica** on the desktop of the station.
- Pour **Lysis Solution MAGNO-sorb**, **Washing Solution 5**, **Washing Solution 6** and **Washing Solution 7** into the appropriate containers placed on the desktop of the station.
- Place the tubes with test samples on the desktop of the station. Choose the rack according to the type of used tubes.
- Add no less than **150 µl of Negative Control (C-)** into the tube for the Negative Control of extraction (C-) and no less than **150 µl of required Positive Control sample** into the tube for the Positive Control of extraction (PCE) (if they are provided for PCR).

NOTE: It is allowed to change the volume and dilution of **Negative Control (C-)** and **Positive Control sample** according to the Instruction manual for the used PCR kit. The volume of control samples should be increased taking into account extra volume.

- Choose and run the protocol of extraction nucleic acid.
- If the elution was carried out into a plate, seal the elution plate with the purified nucleic acid using a self-adhesive film at the end of the extraction if necessary.

Storage of purified nucleic acids

Purified DNA samples can be stored:

- at the temperature from 2 to 8 °C for 1 week;
- at the temperature not more than minus 16 °C for 1 year.

Purified RNA samples can be stored:

- at the temperature from 2 to 8 °C for 4 hours;
- at the temperature from minus 24 to minus 16 °C for 1 month;
- at the temperature not more than minus 68 °C for 1 year.

In case of manual extraction transfer the supernatant for purified nucleic acids storage into a new tube without disturbing the sediment.

In case of automatic extraction transfer eluate into the tubes for purified nucleic acids storage at subzero temperature.

8.3 Variant 100-100M

The volume of test sample for the extraction is 100 µl.

Test material:

- blood plasma,
- nasopharyngeal and oropharyngeal swabs,
- sputum / endotracheal aspirate,
- bronchoalveolar lavage fluid / bronchial washings,
- tissue (biopsy, surgical, autopsy) native material,
- anal canal/rectal mucosal discharge (swab, scrape).

NOTE: See the information about pretreatment of biological material for DNA/RNA extraction in the Instruction manual for the used PCR kit.

NOTE: To work with RNA, it is necessary to use only disposable sterile plastic consumables that have a special marking RNase-free, DNase-free.

NOTE: **Lysis Solution MAGNO-sorb** has unpleasant odor. Work in the PCR box.

Manual procedure when using magnetic racks

- If there is a precipitate in **Lysis Solution MAGNO-sorb** and **Washing Solution 5** warm them up at 60 °C until crystals disappear. Resuspend **Magnetized silica** on vortex.
- Prepare the required number of disposable 1.5-ml tubes (including a tube for the Negative Control of extraction (C) and a tube for the Positive Control of Extraction (PCE) if provided with the amplification kit). Mark the tubes.
- Mix in a disposable tube **Internal Control (IC)** (if it is provided for analysis of given infectious agent), **Component A**, and **Magnetized silica** in the following proportion calculated per one sample: **10 µl of Internal Control (IC)**, **5 µl of Component A**, and **10 µl of Magnetized Silica**. Do not forget to add extra volumes for one more reaction. For example:

| Number of samples for DNA/RNA extraction | Internal Control (IC), µl | Component A, µl | Magnetized silica, µl |
|--|---------------------------|-----------------|-----------------------|
| 6 | 70 | 35 | 70 |
| 12 | 130 | 65 | 130 |
| 18 | 190 | 95 | 190 |
| 24 | 250 | 125 | 250 |

NOTE: It is allowed to change the volume and dilution of **Internal Control (IC)** according to the Instruction manual for the used PCR kit

- Add **25 µl** of the resuspended mixture of **Internal Control (IC)**, **Component A**, and **Magnetized silica** into each tube (recalculate the volume of the mixture if the volume of Internal Control (IC) has been changed);
- Add **300 µl of Lysis Solution MAGNO-sorb** into the tubes.
- Add **100 µl of test sample** into each tube, using a separate filter tip for each sample. Mix by vortexing.

NOTE: It is allowed to change the test sample volumes according to the Instruction manual for the used PCR kit.

- Add **100 µl of the Negative Control (C-)** reagent into the tube for the Negative Control of extraction (C-) and **100 µl of required Positive Control sample** into the tube for the Positive Control of extraction (PCE) (if they are provided for PCR). Tightly close the tubes and thoroughly mix by vortexing.

NOTE: It is allowed to change the volume and dilution of **Negative Control (C-)** and **Positive Control** according to the Instruction manual for the used PCR kit

- Incubate the tubes at **60 °C for 10 min**.
- Centrifuge the tubes shortly, transfer the tubes to a magnetic rack, and incubate for **5 min**. If it is inconvenient/impossible to open the tubes in a magnetic rack without sorbent resuspension, it is necessary to open the tubes before placing them in a magnetic rack.

NOTE: If it is inconvenient/impossible to open the tubes in a magnetic rack without sorbent resuspension, it is necessary to open the tubes before placing them in a magnetic rack.

- Leaving the tubes in the magnetic rack, carefully remove the supernatant inserting a tip near the internal tube wall and using vacuum aspirator. Take a new 200-µl tip for each sample.
- Add **500 µl of Washing Solution 5** to the tubes. Tightly cap the tubes.
- If it is inconvenient/impossible to close tightly the tubes in a magnetic rack, it is necessary to place the tubes in a regular tube rack.
- Resuspend **Magnetized silica** by mixing on vortex. Then sediment the drops by short centrifugation.

13. Transfer the tubes to a magnetic rack for **5 min**.
14. Remove the supernatant as described above.
15. Add **500 µl of Washing Solution 4** to the tubes. Tightly cap the tubes.
16. Resuspend **Magnetized silica** by mixing on vortex. Then sediment the drops by short centrifugation.
17. Transfer the tubes to a magnetic rack for **2 min**.
18. Remove supernatant as described above.
19. Dry the magnetized silica by placing the test tubes with open caps in a thermostat at **60 °C** for **5 minutes**. In case of incomplete drying of the sorbent, the time may be increased up to **10 minutes**.
20. Add **100 µl of Buffer for elution A** to each tube, close the tubes and mix on vortex.

NOTE: It is allowed to change the elution volume according to the Instruction manual for the used PCR kit.

21. Incubate the tubes at **60 °C** for **5 min**. Vortex the tubes in **2 min**.
 22. Sediment the drops by short centrifugation and transfer the tubes to the magnetic rack for **2 min**.
- Supernatant contains purified DNA and RNA.

NOTE: Do not take the tubes away from a magnetic rack when removing RNA/DNA.

Procedure for automated stations for nucleic acid extraction

The detailed information work with automated stations see in the *Instruction Manual* given by the manufacturer of the instrument.

Extraction using stations with magnetic stirring (for example, Auto-Pure 96 (Hangzhou Allsheng Instruments Co., Ltd., China), KingFisher Flex (Thermo FS, Finland)), Nexor 96 (Yantai Addcare Bio-Tech Co., Ltd., China)).

Only the protocols located on the amplisens.ru website should be used for extraction.

NOTE: The extraction protocol for the automated station contains specific parameters that are not available for manual extraction.

1. If there is a precipitate in **Lysis Solution MAGNO-sorb** and **Washing Solution 5** warm them up at **60 °C** until crystals disappear. Resuspend **Magnetized silica** on vortex.
 2. Prepare the required number of 96-well plates (see point 3) compatible with the instrument; mark them according to the protocol.
 3. Add in the plates according to the marking and quantity of samples for extraction:
 - **300 µl of Lysis Solution MAGNO-sorb** in a deepwell plate,
 - **500 µl of Washing Solution 5** in a deepwell plate,
 - **500 µl of Washing Solution 4** in a deepwell plate,
 - **100 µl of Buffer for elution A** in the plate for elution.
- NOTE:** It is allowed to change the elution volume according to the Instruction manual for the used PCR kit.
4. Mix in a separate tube in the following proportion calculated per one sample:
 - **10 µl of Internal Control (IC)** (if it is provided for analysis),
 - **5 µl of Component A**,
 - **10 µl of Magnetized silica**.

Prepare mixture for total number of test and control samples plus extra samples (for example, in case of extraction from 96 samples it is recommended to prepare mixture for 100 extractions, i.e. plus 4 extra samples).

NOTE: It is allowed to change the volume of **Internal Control (IC)** according to the Instruction manual for the used PCR kit.

5. Add **25 µl** of thoroughly resuspended mixture to each well with **Lysis Solution MAGNO-sorb** (recalculate the volume of mixture if the volume of Internal Control (IC) has been changed).
 6. Add **100 µl of test and control** (if they are provided for analysis) **samples** in each well with **Lysis Solution MAGNO-sorb** and mixture according to the marking, using a separate filter tip for each sample.
- NOTE:** It is allowed to change the test and control sample volumes according to the Instruction manual for the used PCR kit.
7. Choose the protocol of extraction, place the prepared plates and tip-plate according to the protocol and placement indicators.
 8. Run the nucleic acid extraction protocol.
 9. If necessary, seal the elution plate with the purified nucleic acid using a self-adhesive film at the end of the extraction.

Storage of purified nucleic acids

Purified DNA samples can be stored:

- at the temperature from 2 to 8 °C for 1 week;
- at the temperature from minus 24 to minus 16 °C for 1 year.

Purified RNA samples can be stored:

- at the temperature from 2 to 8 °C for 4 hours;
- at the temperature from minus 24 to minus 16 °C for 1 month;
- at the temperature not more than minus 68 °C for 1 year.

In case of manual extraction transfer the supernatant for purified nucleic acids storage into a new tube without disturbing the sediment.

In case of automatic extraction transfer eluate into the tubes for purified nucleic acids storage at subzero temperature.

8.4 Variant 100-200M

The volume of test sample for the extraction is 200 µl.

Test material:

- blood plasma,
- whole blood,
- cerebrospinal fluid,
- discharge of the conjunctiva,
- saliva,
- nasopharyngeal and oropharyngeal swabs,
- sputum / endotracheal aspirate,
- bronchoalveolar lavage fluid / bronchial washings,
- pleural/ascitic fluid,
- pus / liquid or dense necrotic content,
- tissue (biopsy, surgical, autopsy) native material,
- paraffin-embedded tissue (biopsy, surgical, autopsy) material,
- urine,
- feces,
- anal canal/rectal mucosal discharge (swab, scrape),
- blood-sucking arthropods (ticks),
- cultures of microorganisms,
- water (drinking water, open sources, sewage),
- washes from environmental objects.

NOTE: See the information about pretreatment of biological material for RNA / DNA extraction in the Instruction manual for the used PCR kit.

NOTE: To work with RNA, it is necessary to use only disposable sterile plastic consumables that have a special marking RNase-free, DNase-free.

NOTE: **Lysis Solution MAGNO-sorb** has unpleasant odor. Work in the PCR box.

Manual procedure when using magnetic racks

1. If there is a precipitate in **Lysis Solution MAGNO-sorb** and **Washing Solution 5** warm them up at **60 °C** until crystals disappear. Resuspend **Magnetized silica** on vortex.
2. Prepare the required number of disposable 1.5-ml tubes (including a tube for the Negative Control of extraction (C) and a tube for the Positive Control of Extraction (PCE) if provided

with the amplification kit). Mark the tubes.

3. Mix in a disposable tube **Internal Control (IC)** (if it is provided for analysis of given infectious agent), **Component A**, and **Magnetized silica** in the following proportion calculated per one sample: **10 µl of Internal Control (IC)**, **10 µl of Component A**, and **20 µl of Magnetized Silica**. Do not forget to add extra volumes for one more reaction. For example:

| Number of samples for DNA/RNA extraction | Internal Control (IC), µl | Component A, µl | Magnetized silica, µl |
|--|---------------------------|-----------------|-----------------------|
| 6 | 70 | 70 | 140 |
| 12 | 130 | 130 | 260 |
| 18 | 190 | 190 | 380 |
| 24 | 250 | 250 | 500 |

NOTE: It is allowed to change the volume and dilution of **Internal Control (IC)** according to the Instruction manual for the used PCR kit.

4. Add **40 µl** of the resuspended mixture of **Internal Control (IC)**, **Component A**, and **Magnetized silica** into each tube (recalculate the volume of mixture if the volume of Internal Control (IC) has been changed);
5. Add **600 µl of Lysis Solution MAGNO-sorb** into the tubes.
6. Add **200 µl of the test sample** into each tube, using a separate filter tip for each sample. Mix by vortexing.

NOTE: It is allowed to change the test sample volumes according to the Instruction manual for the used PCR kit.

7. Add **200 µl of the Negative Control (C-)** reagent into the tube for the Negative Control of extraction (C-) and **200 µl of required Positive Control sample** into the tube for the Positive Control of extraction (PCE) (if they are provided for PCR). Tightly close the tubes and thoroughly mix by vortexing.

NOTE: It is allowed to change the volume and dilution of **Negative Control (C-)** and **Positive Control** according to the Instruction manual for the used PCR kit.

8. Incubate the tubes at **60 °C** for **10 min**.
9. Centrifuge the tubes shortly, transfer the tubes to a magnetic rack, and incubate for **5 min**. If it is inconvenient/impossible to open the tubes in a magnetic rack without sorbent resuspension, it is necessary to open the tubes before placing them in a magnetic rack.

NOTE: Leaving the tubes in the magnetic rack, carefully remove the supernatant inserting a tip near the internal tube wall and using vacuum aspirator. Take a new 200-µl tip for each sample.

11. Add **500 µl of Washing Solution 5** to the tubes. Tightly cap the tubes.

NOTE: If it is inconvenient/impossible to close tightly the tubes in a magnetic rack it is necessary to place the tubes in a regular tube rack.

12. Resuspend **Magnetized silica** by mixing on vortex. Then sediment the drops by short centrifugation.
 13. Transfer the tubes to a magnetic rack for **5 min**.
 14. Remove the supernatant as described above.
 15. Add **500 µl of Washing Solution 4** to the tubes. Tightly cap the tubes.
 16. Resuspend **Magnetized silica** by mixing on vortex. Then sediment the drops by short centrifugation.
 17. Transfer the tubes to a magnetic rack for **2 min**.
 18. Remove supernatant as described above
 19. Dry the **Magnetized silica** by placing the test tubes with open caps in a thermostat at **60 °C** for **5 minutes**. In case of incomplete drying of the sorbent, the time may be increased up to **10 minutes**.
 20. Add **100 µl of Buffer for elution A** to each tube, close the tubes and mix on vortex.
- NOTE:** It is allowed to change the elution volume according to the Instruction manual for the used PCR kit.
21. Incubate the tubes at **60 °C** for **5 min**. Vortex the tubes in **2 min**.
 22. Sediment the drops by short centrifugation and transfer the tubes to the magnetic rack for **2 min**.
- Supernatant contains purified RNA and DNA.

NOTE: Do not take the tubes away from a magnetic rack when removing RNA/DNA.

Procedure for automated stations for nucleic acid extraction

The detailed information work with automated stations see in the *Instruction Manual* given by the manufacturer of the instrument.

Extraction using stations with magnetic stirring (for example, Auto-Pure 96 (Hangzhou Allsheng Instruments Co., Ltd., China), KingFisher Flex (Thermo Fisher Scientific Oy, Finland), Nexor 96 (Yantai Addcare Bio-Tech Co., Ltd., China)).

Only the protocols located on the amplisens.ru website should be used for extraction.

NOTE: The extraction protocol for the automated station contains specific parameters that are not available for manual extraction

1. If there is a precipitate in **Lysis Solution MAGNO-sorb** and **Washing Solution 5** warm them up at **60 °C** until crystals disappear. Resuspend **Magnetized silica** on vortex.
 2. Prepare the required number of 96-well plates (see point 3) compatible with the instrument; mark them according to the protocol.
 3. Add in the plates according to the marking and quantity of samples for extraction:
 - **600 µl of Lysis Solution MAGNO-sorb** in a deepwell plate,
 - **500 µl of Washing Solution 5** in a deepwell plate,
 - **500 µl of Washing Solution 4** in a deepwell plate,
 - **100 µl of Buffer for elution A** in the plate for elution.
- NOTE:** It is allowed to change the elution volume according to the Instruction manual for the used PCR kit.
4. Mix in a separate tube in the following proportion calculated per one sample:
 - **10 µl of Internal Control (IC)** (if it is provided for analysis),
 - **10 µl of Component A**,
 - **20 µl of Magnetized silica**.

Prepare mixture for total number of test and control samples plus extra samples (for example, in case of extraction from 96 samples it is recommended to prepare mixture for 100 extractions, i.e. plus 4 extra samples).

NOTE: It is allowed to change the volume and dilution of **Internal Control (IC)** according to the Instruction manual for the used PCR kit.

5. Add **40 µl** of thoroughly resuspended mixture to each well with **Lysis Solution MAGNO-sorb** (recalculate the volume of mixture if the volume of Internal Control (IC) has been changed).
 6. Add **200 µl of test and control** (if they are provided for analysis) **samples** in each well with **Lysis Solution MAGNO-sorb** and mixture according to the marking, using a separate filter tip for each sample.
- NOTE:** It is allowed to change test and control sample volumes according to the Instruction manual for the used PCR kit.
7. Choose the protocol of extraction, place the prepared plates and tip-plate according to the protocol and placement indicators.
 8. Run the nucleic acid extraction protocol.
 9. If necessary, seal the elution plate with the purified nucleic acid using a self-adhesive film at the end of the extraction.

Storage of purified nucleic acids

Purified DNA samples can be stored:

- at the temperature from 2 to 8 °C for 1 week;
- at the temperature from minus 24 to minus 16 °C for 1 year.

Purified RNA samples can be stored:

- at the temperature from 2 to 8 °C for 4 hours;
- at the temperature from minus 24 to minus 16 °C for 1 month;
- at the temperature not more than minus 68 °C for 1 year.

In case of manual extraction transfer the supernatant for purified nucleic acids storage into a new tube without disturbing the sediment.

In case of automatic extraction transfer eluate into the tubes for purified nucleic acids storage at subzero temperature.

9. TROUBLESHOOTING

These troubleshooting rules may be helpful in explaining any questions that may arise.

False negatives with extraction product:

- Degradation of the nucleic acid contained in the sample. It's necessary to use a new sample, store samples appropriately.
- Loss of nucleic acid pellet. Carefully draw off the washing solution and try not to remove the sorbent.
- Degradation of the extracted nucleic acid. It's necessary to use plastic free from DNases and RNases.

False positives with extraction product:

- Contamination during sample extraction. It's necessary to open one test tube at a time. Avoid spilling the contents of the test tube, always change tips.
- Contamination of the reagents prepared for the step. It's necessary to repeat the test.
- Contamination of the Extraction Area by amplicons. It's necessary to clean surfaces and instruments using aqueous detergents, wash lab coats, replace test tubes and tips in use. Use different laboratory coats in Extraction, Amplification and Detection areas.

If you have any further questions or if encounter problems, please contact our Authorized representative in the European Community.

10. TRANSPORTATION

MAGNO-sorb nucleic acid extraction kit should be transported at 2–25 °C for no longer than 5 days.

11. STABILITY AND STORAGE

All components of **MAGNO-sorb** nucleic acid extraction kit **variant 100-200, variant 100-100M, variant 100-200M and variant 100-1000** are to be stored at 2–25 °C, when not in use. All components of **MAGNO-sorb** nucleic acid extraction kit are stable until the expiry date stated on the label. The shelf life of reagents before and after the first use is the same, unless otherwise stated.

12. REFERENCES

1. Karpov S.I., Matveeva M.V., Selemenev V.F. // Russian Journal of Physical Chemistry A. 2001. Vol. 75, No 2, pp. 266-271;
2. Sameer A., Barghouthi A // Indian J. Microbiol. 2011. Vol. 51, pp. 430-444;
3. Chiang C.-L., Sung C.-S., Wu T.-F., et al. Application of Superparamagnetic Nanoparticles in Purification of Plasmid DNA from Bacterial Cells // J. Chromatogr. B. 2005. V. 822. P. 54-60.

13. QUALITY CONTROL

In compliance with Federal Budget Institute of Science "Central Research Institute for Epidemiology" ISO 13485-Certified Quality Management System, each lot of the **MAGNO-sorb** nucleic acid extraction kit has been tested against predetermined specifications to ensure consistent product quality.

Please contact our Authorized representative in the European Community if side effects, undesirable reactions, facts and circumstances that pose a threat to the life and health of citizens and medical workers are identified during the use of the reagent kit.

List of Changes Made in the Instruction Manual

| VER | Location of changes | Essence of changes |
|---------------------------|---|---|
| 29.09.17 ME | 5. General precautions, 14. Key to symbols used | Information about hazards was rewritten according to the Regulation 1272/2008/EC |
| 26.02.19 PM | 3. Content | The colour of the reagents was specified |
| 08.04.20 KK | Through the text | The text formatting was changed |
| | Footer | The phrase "Not for use in the Russian Federation" was added |
| 11.03.21 VA | — | The name, address and contact information for Authorized representative in the European Community was changed |
| 08.06.22 MM | 1. Intended use | The intended use was specified. The list of biological material was expanded. The "Indications and contraindications for use of the reagent kit" subsection was added |
| | Through the text | Component A was changed to Component A-2 for variant 100-1000. Buffer for elution B was changed to Buffer for elution for variant 100-1000 and for variant 100-200. The variant 100-100M and variant 100-200M were added |
| | 3. Content | Volumes and quantity of reagents were changed for variant 100-200 and for variant 100-1000. REF K3-1061-100-CE and REF K3-1062-100-CE were added for variant 100-200 and variant 100-200 respectively |
| | 4. Additional requirements | The automatic stations for nucleic acid extraction and additional materials for their use were added |
| | 5. General precautions | Hazard identification information for Lysis Solution MS-URO, Washing solution 4 and Mucolysin reagents were added in accordance with Regulation 1272/2008/EC |
| | 6. Sampling and handling | "Interfering substances and limitations of using test material samples" and "Potential interfering substances" subsections were added |
| | 8. Protocol | The procedure of manual extraction using magnetic racks was clarified for variant 100-200 and for variant 100-1000. The procedure of nucleic acid extraction using automated stations was added for variant 100-200 and for variant 100-1000. The procedure of nucleic acid extraction were added for variant 100-100M and variant 100-200M |
| | 10. Transportation | The storage temperature of transportation was changed from 2-8 °C to 2-25 °C |
| | 13. Quality control | The Authorized representative in the European Community was specified for the contact in case of undesirable effects when using the reagent kit |
| | 08.08.23 BA | 3. Content |
| Footer | | The reference numbers K2-16-200-CE, K2-16-1000-CE were deleted |
| 30.01.25 HM | 1. Intended use | The intended use was specified. The list of biological material was reworded |
| | 3. Content | For variant 100-1000 Component A-2 was changed to the reagent Component A. For all forms Buffer for elution was renamed to Buffer for elution A |
| | 4. Additional requirements | Auto-Pure 24, Nexor 96 automated DNA extraction systems were added. NEON 100 DNA extraction pipetting system was renamed to Xiril (Neon 100 series) |
| | 5. General precautions | Hazard identification information for Washing Solution 4 and Washing Solution 5 was specified |
| | 6. Sampling and handling | The list of biological material was specified |
| | 8. Protocol | Additional clarifying information was added |
| 11. Stability and storage | Storage temperature for variant 100-1000 was changed from 2–8 °C to 2–25 °C | |

AmpliSens®



Ecoli Dx, s.r.o., Purkyňova 74/2
110 00 Praha 1, Czech Republic
Tel.: +420 325 209 912
Cell: +420 739 802 523



Federal Budget Institute of
Science "Central Research
Institute for Epidemiology"
3A Novogireevskaya Street
Moscow 111123 Russia