

AmpliSens® HBV-Monitor-FRT PCR kit



For Professional Use Only

Instruction Manual

KEY TO SYMBOLS USED

	Catalogue number		Caution
	Batch code		Contains sufficient for <n> tests
	Research use only		Use-by Date
	Version		Consult instructions for use
	Temperature limit		Keep away from sunlight
	Manufacturer	NCA	Negative control of amplification
	Date of manufacture	C-	Negative control of extraction
PCE	Positive Control of Extraction	C+	Positive control of Amplification
		IC	Internal control

1. INTENDED USE

AmpliSens® HBV-Monitor-FRT PCR kit is not a medical device. PCR kit is an *in vitro* nucleic acid amplification test for quantitative detection of hepatitis B virus (HBV) DNA in biological material (blood plasma) using real-time hybridization-fluorescence detection of amplified products.

NOTE: For research use only. Not for diagnostic procedures.

2. PRINCIPLE OF PCR DETECTION

Principle of testing is based on the DNA extraction from the samples of test material with the exogenous internal control sample (Internal Control STI-87 (IC)) and simultaneous amplification of DNA fragments with hybridization-fluorescence detection. Exogenous internal control (Internal Control STI-87 (IC)) allows to control all PCR-analysis stages of each individual sample and to identify possible reaction inhibition.

AmpliSens® HBV-Monitor-FRT PCR kit uses "hot-start", which greatly reduces frequency of nonspecifically primed reactions. "Hot-start" is guaranteed by using chemically modified polymerase (TaqF). Chemically modified polymerase (TaqF) is activated by heating at 95 °C for 15 min.

The PCR kit contains the system for prevention of contamination by amplicons using the enzyme uracil-DNA-glycosylase (UDG) and deoxyuridine triphosphate (dUTP). The results of amplification are registered in the following fluorescence channels:

Table 1

Channel for fluorophore	FAM	JOE
DNA-target	IC DNA	HBV DNA
Target gene	Artificially synthesized sequence	Precore/core

3. CONTENT

AmpliSens® HBV-Monitor-FRT PCR kit is produced in 2 forms:

Form 1: PCR kit variant FRT, Negative Control (C-) (4 tubes), HBV-Q calibration kit, R-V5-MC(RG,iQ,Mx,Dt)-CE,

Form 2: PCR kit variant FRT, Negative Control (C-) (4 tubes), HBV-Q calibration kit in bulk¹, R-V5-MC(RG,iQ,Mx,Dt)-CE-B.

PCR kit variant FRT includes:

Reagent	Description	Volume, ml	Quantity	
PCR-mix-1-FL HBV	clear liquid from colorless to light lilac colour	0.3	4 tubes	
PCR-mix-2-FRT	colorless clear liquid	0.2	4 tubes	
Polymerase (TaqF)	colorless clear liquid	0.02	4 tubes	
DNA calibrators	PIC1 HBV	colorless clear liquid	0.1	4 tubes
	PIC2 HBV	colorless clear liquid	0.1	4 tubes
Buffer for elution	colorless clear liquid	1.2	4 tubes	
Positive Control-1-HBV*	colorless clear liquid	0.06	4 tubes	
Positive Control-2-HBV*	colorless clear liquid	0.06	4 tubes	
Internal Control STI-87 (IC)*	colorless clear liquid	0.28	4 tubes	

* must be used in the extraction procedure (see section 8.1 for details).

PCR kit variant FRT is intended for 80 reactions, including controls and DNA calibrators.

The following reagent is additionally included in the PCR kit variant FRT:

Reagent	Description	Volume, ml	Quantity
Negative Control (C-)*	colorless clear liquid	1.2	4 tubes

* must be used in the extraction procedure as Negative Control of Extraction.

HBV-Q calibration kit includes:

Reagent	Description	Volume, ml	Quantity
Calibrator HBV-Q	yellow powder	-	1 tube
Solvent Q	colourless clear liquid	1.2	3 tubes

PCR kit also includes:

The software in Microsoft® Excel format for processing of data and generation of results.

4. ADDITIONAL REQUIREMENTS

- DNA/RNA extraction kit.
- Sterile RNase-free/DNase-free pipette tips with filters up to 200 µl and up to 1000 µl.
- Tube racks.
- Vortex mixer.
- Desktop microcentrifuge up to 12,000 g (suitable for Eppendorf tubes).
- PCR box.
- Real-time instruments (for example, Rotor-Gene 3000/6000 (Corbett Research, Australia); iCycler iQ5 (Bio-Rad, USA); Mx3000P (Stratagene, USA)).
- Disposable polypropylene tubes:
 - a) thin-walled 0.2-ml PCR tubes with optical transparent domed or flat caps or strips of eight 0.2-ml tubes with optical transparent caps if a plate-type instrument is used;
 - b) thin-walled 0.2-ml PCR tubes with flat caps or strips of four 0.1-ml Rotor-Gene PCR tubes if a rotor-type instrument is used.
- Pipettes (adjustable).
- Refrigerator for 2-8 °C.
- Deep-freezer at the temperature from minus 24 to minus 16 °C.
- Reservoir for used tips.
- Disposable powder-free gloves and a laboratory coat.

5. GENERAL PRECAUTIONS

The user should always pay attention to the following:

- Use sterile pipette tips with aerosol filters and use a new tip for every procedure.
- Store all extracted positive material (specimens, controls and amplicons) away from all other reagents and add it to the reaction mix in a distantly separated facility.
- Thaw all components thoroughly at room temperature before starting an assay.
- When thawed, mix the components and centrifuge briefly.
- Use disposable protective gloves and laboratory cloths, and protect eyes while samples and reagents handling. Thoroughly wash hands afterwards.
- Do not eat, drink, smoke, apply cosmetics, or handle contact lenses in laboratory work areas.
- Do not use the PCR kit if the internal packaging was damaged or its appearance was changed.
- Do not use the PCR kit if the transportation and storage conditions according to the Instruction Manual were not observed.
- Do not use a kit after its expiration date.
- Dispose of all specimens and unused reagents in accordance with local regulations.
- Samples should be considered potentially infectious and handled in biological cabinet in compliance with appropriate biosafety practices.
- Clean and disinfect all samples or reagents spills using a disinfectant, such as 0.5 % sodium hypochlorite or another suitable disinfectant.
- Avoid inhalation of vapors, samples and reagents contact with the skin, eyes, and mucous membranes. Harmful if swallowed. If these solutions come into contact, rinse the injured area immediately with water and seek medical advice if necessary.
- While observing the conditions of transportation, operation and storage, there are no risks of explosion and ignition.
- Safety Data Sheets (SDS) are available on request.
- The PCR kit is intended for single use for PCR analysis of specified number of samples (see the section "Content").
- The PCR kit is ready for use in accordance with the Instruction Manual. Use the PCR kit strictly for intended purpose.
- Use of this product should be limited to personnel trained in DNA amplification techniques.
- Workflow in the laboratory must be one-directional, beginning in the Extraction Area and moving to the Amplification and Detection Area. Do not return samples, equipment and reagents in the area where the previous step was performed.



Some components of this kit contain sodium azide as a preservative. Do not use metal tubing for reagent transfer.

¹ In bulk form contains unlabeled tubes. Tubes with identical reagent are packed in one bag with label.

6. SAMPLING AND HANDLING

NOTE: Obtaining biological material samples for PCR-analysis, transportation, and storage are described in manufacturer's handbook [1]. It is recommended that this handbook is read before starting work.

AmpliSens® HBV-Monitor-FRT PCR kit is intended for the analysis of the DNA extracted with DNA extraction kits from the biological material (blood plasma):
Blood samples are collected in the morning on an empty stomach into the tube with EDTA solution as the anticoagulant. Several times invert the closed tubes to ensure proper mixing. To collect plasma, centrifuge the tubes with blood at 800–1600 g for 20 min for 6 h after blood taking. Remove obtained plasma and transfer to the new tubes.
The blood serum may also be used in some cases. The analytical sensitivity of the reagent kit is retained for this material; however, the clinical sensitivity may be significantly decreased as a result of viral particles precipitation during blood clot retraction.

NOTE: The peripheral blood plasma and the blood serum samples can be stored for no longer than 3 days at a temperature range from 2 to 8 °C or, if stored at a temperature no greater than minus 68 °C - for a long time.

7. WORKING CONDITIONS

AmpliSens® HBV-Monitor-FRT PCR kit should be used at the temperature from 20 to 28 °C and relative humidity from 15 to 75 %.

8. PROTOCOL

8.1 DNA extraction

It is recommended to use the following nucleic acid extraction kits:

- **RIBO-prep**, **REF** K2-9-Et-100-CE
- **MAGNO-sorb**, **REF** K2-16-1000-CE;
- NucliSENS easyMAG automated nucleic acid extraction system (bioMérieux, France) can also be used (see Guidelines [2] for details).

If using RIBO-prep kit extract the RNA/DNA according to the manufacturer's protocol taking into account next additions and improvements

- If a large number of samples is being tested, it is acceptable to mix the **Solution for Lysis** and the **Internal Control** in a separate sterile flask (based on addition of **300 µl of Solution for Lysis** and **10 µl of Internal Control** per one sample), followed by a transfer of **300 µl** of the prepared mix into each of the previously prepared **1.5 ml** tubes.

NOTE: For each panel it is necessary to set up **two** positive controls of extraction – **PCE-1** and **PCE-2**. To the tube labelled **PCE-1** add **90 µl** of **Negative Control** and **10 µl** of **Positive Control-1-HBV**, to the tube labelled **PCE-2** add **90 µl** of **Negative Control** and **10 µl** of **Positive Control-2-HBV**. Close the lids and vortex the tubes.

• Centrifuge at 12,000 g throughout the extraction procedure.

NOTE: **If using the MAGNO-sorb kit** extract the RNA/DNA according to the manufacturer's protocol taking into account following additions and improvements:

- In case of DNA extraction from blood plasma sample of 1000 µl, the volume of the **Internal Control STI-87** required for 24-tube panel is **0.28 ml**. In case of other panels and DNA extraction from blood plasma sample of 200 µl see the **MAGNO-sorb** instruction manual.
- To prepare the Positive Control of Extraction 1 (**PCE-1**), add **90 µl** of the **Negative Control (C-)** sample and **10 µl** of the **Positive Control-1-HBV** sample to the new tube containing **Lysis Solution MAGNO-sorb**.
- To prepare the Positive Control of Extraction 2 (**PCE-2**), add **90 µl** of the **Negative Control (C-)** sample and **10 µl** of the **Positive Control-2-HBV** sample to the new tube containing **Lysis Solution MAGNO-sorb**.
- To prepare the Negative Control of Extraction (**C-**), add **100 µl** of the **Negative Control (C-)** sample to the new tube containing **Lysis Solution MAGNO-sorb**.
- The volume of **Buffer for elution** required for extraction from both 1000 and 200 µl of blood plasma samples is **70 µl**.

8.1.1 Calibration and calculation of the coefficient B using **HBV-Q** calibration kit

If coefficient B for the extraction kit/automatic platform is not specified in the *Important Product Information Bulletin* enclosed in the PCR kit, the calibration for calculation of coefficient B should be carried out by oneself with the aid of the **HBV-Q** calibration kit included in this PCR kit. See below for details.

NOTE:

The calibration procedure is necessary to define Coefficient B and it is performed during the **first PCR run** for the given lot. Calibration is performed **only once** for each new lot of the **AmpliSens® HBV-Monitor-FRT** PCR kit and is conducted using the same DNA extraction kit/automatic station as was used in the PCR assay. To carry out calibration, it is necessary to analyse 5 extra samples: the repeat of **Positive Control-1-HBV**, the repeat of **Positive Control-2-HBV**, and calibrator **HBV-Q** in triplicate.

Calibrator **HBV-Q** preparation

1. Vortex the tube with calibrator **HBV-Q**, gently open the tube, and add **400 µl** of **solvent Q** avoiding the contents spraying. Use tips with filters.
2. Close the tube and incubate it at room temperature for 20 min vortexing periodically.
3. Once the contents are fully dissolved, vortex the tube for 3-5 s to make sure that there are no drops on the caps of the tube.

Perform calibration with the same DNA extraction kit used in the PCR assay.

NOTE: Extract the DNA according to the manufacturer's protocol.

Transfer 10 µl of **Internal Control STI-87 (IC)** (per one sample) to samples or **Lysis solution** before extraction.

In case of extracting from **100 µl of plasma**, add dissolved calibrator **HBV-Q** to three tubes for DNA extraction (100 µl per each tube).

In case of extracting from any **other plasma volume (100 – 1000 µl)**, transfer dissolved calibrator **HBV-Q** to three tubes for DNA extraction (100 µl per each tube) and add **solvent Q** up to the extraction volume (for example, if the extraction volume is 1 ml then add 100 µl of calibrator **HBV-Q** and 900 µl of solvent Q).

When extraction is completed, perform PCR as described in this *Instruction manual*.

Use the mean concentration values obtained in the channels for the **FAM** and **JOE** fluorophores for three repeats with calibrator **HBV-Q** for calculation of coefficient B using the following formula:

$$\text{Coefficient B} = \frac{\text{IC DNA copies in calibrator HBV-Q (FAM channel)}}{\text{HBV DNA copies in calibrator HBV-Q (JOE channel)}} \times \text{coefficient C}$$

Coefficient C is specified in the *Important Product Information Bulletin* enclosed in the **AmpliSens® HBV-Monitor-FRT** PCR kit.

NOTE: The calculated value of coefficient B should be within range specified in the *Important Product Information Bulletin* enclosed in the applied PCR kit lot
Write down the coefficient B value in the *Important Product Information Bulletin* enclosed with the given lot of the PCR kit and use it for concentrations calculation of biological and control samples (See Data Analysis section).

Also see the Guidelines [2] to **AmpliSens® HBV-Monitor-FRT** PCR kit.

Write down the calculated values for **Positive Control-1-HBV** and **Positive Control-2-HBV** in the *Important Product Information Bulletin* enclosed with the given lot of the PCR kit. Determine the mean value for both **Positive Control-1-HBV** and for **Positive Control-2-HBV**. Set the acceptable value range for both **Positive Control-1-HBV** and for **Positive Control-2-HBV** as follows: from "calculated mean value"/3 to "calculated mean value"x3.

For example,
the calculated values for **Positive Control-1-HBV** in two replicates are 500,000 IU/ml and 700,000 IU/ml;
the calculated mean value for **Positive Control-1-HBV** is 600,000 IU/ml;
the acceptable value range for **Positive Control-1-HBV** varies from 200,000 to 1800,000 IU/ml.

Write down the calculated acceptable value range for **Positive Control-1-HBV** and for **Positive Control-2-HBV** in the *Important Product Information Bulletin*, and use it to verify further assays conducted using this lot of the PCR kit (see 9. Data Analysis section).

8.2 Preparing PCR

The total reaction volume is **50 µl**, the volume of DNA sample is **25 µl**.

8.2.1 Preparing tubes for PCR

NOTE: All components of the reaction mixture should be mixed just before use. See Table 2 for the reaction mixture preparation scheme.

1. Before starting work, thaw and thoroughly vortex all reagents of the kit. Make sure that there are no drops on the caps of the tubes.
2. Take the required number of the PCR tubes for amplification for the biological and control samples (including 3 controls of extraction and 4 calibrators). The type of tubes depends on the PCR instrument used for analysis.
3. To prepare the reaction mixture, mix reagents **15 µl of PCR-mix-1-FL HBV**, **10 µl of PCR-mix-2-FRT** and **1.0 µl of polymerase (TaqF)** per one reaction in a new sterile tube. Thoroughly vortex the mixture, make sure that there are no drops on the caps of the tubes.

It is recommended that the reaction mixture for 20 reactions is prepared in case of extraction from 12 to 16 samples (two NucliSENS easyMAG arrays). To do this, into the tube with **PCR-mix-1-FL HBV** transfer the entire contents of the tube with **PCR-mix-2-FRT** and the entire contents of the tube with **polymerase (TaqF)**. **Do not store the prepared mixture!**

Table 2
Scheme of reaction mixture preparation

Reagent volume per one reaction, µl			Reagent volumes per number of samples specified plus 1 extra reaction		
			15.00	10.00	1.00
Number of biological samples	Number of extracted samples ²	Number of samples analysed in PCR ³	PCR-mix-1-FL HBV	PCR-mix-2-FRT	Polymerase (TaqF)
3	6	10	165	110	11
4	7	11	180	120	12
5	8 ⁴	12	195	130	13
6	9	13	210	140	14
7	10	14	225	150	15
8	11	15	240	160	16
9	12 ⁵	16	255	170	17
10	13	17	270	180	18
11	14	18	285	190	19
12	15	19	Entire contents of the tube	Entire contents of the tube	Entire contents of the tube
13	16 ⁶	20	Entire contents of the tube	Entire contents of the tube	Entire contents of the tube

4. Transfer **25 µl** of the prepared reaction mixture to each PCR tube. Discard the unused reaction mixture.
5. Using tips with filters add **25 µl** of biological **DNA-samples**. Thoroughly mix by pipetting. Avoid formation of air bubbles.

NOTE: Avoid transferring of the sorbent together with the DNA sample in case of extraction with NucliSENS easyMAG automated nucleic acid extraction system or **MAGNO-sorb** kit.

6. Carry out the control amplification reactions:
 - PCE-1** – Add **25 µl** of **DNA sample extracted from Positive Control-1-HBV** to the tube labelled **PCE-1** (Positive Control of Extraction).
 - PCE-2** – Add **25 µl** of **DNA sample extracted from Positive Control-2-HBV** to the tube labelled **PCE-2** (Positive Control of Extraction).
 - C-** – Add **25 µl** of **DNA sample extracted from Negative Control (C-)** to the tube labelled **C-** (Negative control of Extraction).
 - C+1** – Add **DNA calibrator PIC1 HBV** to the two tubes labelled **C+1** (Positive Control of Amplification) (**25 µl per each tube**).
 - C+2** – Add **DNA calibrator PIC2 HBV** to the two tubes labelled **C+2** (Positive Control of Amplification) (**25 µl per each tube**).
- NOTE:** It is also necessary to carry out **Negative Control of Amplification (NCA)** at suspicion on possible contamination.
- NCA** – Add **25 µl** of **Buffer for elution** to the tube labelled **NCA** (Negative Control of Amplification).

NOTE: Mix the tubes thoroughly by pipetting avoiding foaming

² Number of biological samples + 3 controls of DNA extraction (N+3, N – number of biological samples).

³ Number of biological samples + 3 controls of DNA extraction + 4 DNA calibrators (N+7, N – number of biological samples).

⁴ Extraction of one strip with NucliSENS easyMAG automated system (8 tubes).

⁵ 12-tube extraction panel.

⁶ Extraction of 2 stripes with NucliSENS easyMAG automated system (16 tubes).

8.2.2 Amplification

1. Create a temperature profile on your instrument as follows:

Table 3

AmpliSens-2 RG program (for rotor-type instruments)				
Step	Temperature, °C	Time	Fluorescence detection	Cycles
1	50	15 min	–	1
2	95	15 min	–	1
3	95	5 s	–	5
	60	20 s	–	
	72	15 s	–	
4	95	5 s	–	40
	60	20 s	FAM, JOE, ROX, Cy5	
	72	15 s	–	

Table 4

AmpliSens-2 iQ program (for plate-type instruments)				
Step	Temperature, °C	Time	Fluorescence detection	Cycles
1	50	15 min	–	1
2	95	15 min	–	1
3	95	5 s	–	5
	60	20 s	–	
	72	15 s	–	
4	95	5 s	–	40
	60	30 s	FAM, HEX, ROX, Cy5	
	72	15 s	–	

NOTE: The use of either AmpliSens-2 RG or AmpliSens-2 iQ programs allows to simultaneously carry out any test combinations just in one run. In case of one instrument simultaneously performing only the DNA *HBV* detection tests, it is possible to omit the first step of the program (50 °C, 15 min) to spare time.

NOTE: Channels for the ROX and Cy5 fluorophores are activated upon request, when the multiplex-format tests are carried out, which use these channels.

- Adjust the fluorescence channel sensitivity according to the *Important Product Information Bulletin* and Guidelines [2].
- Insert the tubes into the reaction module of the device.
- Run the amplification program with fluorescence detection.
- Analyse results after the amplification program is completed.

9. DATA ANALYSIS

Analysis of results is performed by the software of the real-time PCR instrument used by measuring the fluorescence signal accumulation in two channels:

- The signal of the Internal Control DNA amplification product is detected in the channel for the FAM fluorophore.
- The signal of the *HBV* DNA fragment amplification product is detected in the channel for the JOE fluorophore.

Results are interpreted by the crossing (or not-crossing) the fluorescence curve with the threshold line set at the specific level (in the middle of the linear fragment of the positive control fluorescence growth in the log scale) that corresponds to the presence (or absence) of the *Ct* values for this sample in the corresponding column of the results grid.

Based on the *Ct* values (the crossing of the fluorescence curve with the threshold line set at the specific level) and on the specified values for the calibrators, PIC1 *HBV* and PIC2 *HBV*, the calibration line is automatically plotted and the values for the number of *HBV* DNA copies (channel for the JOE fluorophore) and for the number of Internal Control DNA copies (channel for the FAM fluorophore) in a PCR sample are calculated. The obtained values are used for the *HBV* DNA concentration calculation in test and control samples, using the formula:

$$\frac{\text{HBV DNA copies per PCR sample}}{\text{IC DNA copies per PCR sample}} \times \text{coefficient A} \times \text{coefficient B} = \text{IU of HBV DNA/ml of plasma}$$

$$\text{Coefficient A} = \frac{100}{\text{extraction volume, } \mu\text{l}}$$

NOTE: Coefficient A = 1 when calculating PCE-1 and PCE-2 concentrations

Coefficient B is specified in the *Important Product Information Bulletin* enclosed to the PCR kit and is specific for each lot. It cannot be used with PCR kits of different lots. Coefficient B should be calculated as the result of calibration during the first PCR run (see section 8.1.1 for details).

If the result is greater than 100,000,000 IU/ml, then it is interpreted as the **greater than 100,000,000 IU/ml result**. If the obtained value is higher than the linear range, then the sample may be re-tested after 10x dilution; the produced result is multiplied by 10.

NOTE: If the result is less than 150 IU/ml (extraction from 100 µl), or less than 75 IU/ml (extraction from 200 µl), or less than 15 IU/ml (extraction from 1 ml), then it is interpreted as the **less than 150, or less than 75, or less than 15 IU/ml result**, respectively.

If it is necessary to obtain results expressed in copies/ml, the results measured in International Units (IU/ml) should be multiplied by 1.7 (1 IU = 1.70 copies, 1 copy = 0.59 IU).

The result of the analysis is considered reliable only if the results obtained for the Positive and Negative Controls of amplification as well as for the Negative and Positive Controls of extraction are correct (see Table 5):

Table 5

Results for controls			
Control	Stage for control	Result of amplification in the channel for fluorophore	
		FAM	JOE
C-	DNA extraction, PCR	Positive (number of IC copies in PCR sample is greater than the boundary value)	Negative
PCE-1	DNA extraction, PCR	Positive (number of IC copies in PCR sample is greater than the boundary value)	Positive (should be within range specified in the Bulletin as a result of calculation with IC copies/ml)
PCE-2	DNA extraction, PCR	Positive (number of IC copies in PCR sample is greater than the boundary value)	Positive (should be within range specified in the Bulletin as a result of calculation with IC copies/ml)
C+1	PCR	Positive	Positive
C+2	PCR	Positive	Positive
NCA	PCR	Negative	Negative

NOTE: Boundary values and the range of values for PCE-1 (Positive Control-1-*HBV*) and PCE-2 (Positive Control-2-*HBV*) calculated with IC copies/ml are specified in the *Important Product Information Bulletin* enclosed to the PCR kit.

10. TROUBLESHOOTING

The results are not taken into account in the following cases:

- The *Ct* value is determined for the Negative Control of Extraction (C-) in the channel for the JOE fluorophore and/or Negative Control of Amplification (NCA) in the channels for the FAM and JOE fluorophores. The PCR analysis (beginning with the DNA extraction stage) should be repeated for all the samples in which *HBV* DNA was detected.
- The concentration of Internal Control in the corresponding channel is less than the boundary value specified in the *Important Product Information Bulletin*. The PCR analysis (beginning with the first analysis stage) should be repeated for the sample.
- The correlation coefficient, R^2 is less than 0.98 when plotting the calibration line, then the amplification and detection for all the samples should be repeated.
- The calculated concentrations of Positive Control-1-*HBV* and Positive Control-2-*HBV* exceed the range specified in the *Important Product Information Bulletin*. The PCR analysis (beginning with the DNA extraction stage) should be repeated for all the samples.

11. TRANSPORTATION

AmpliSens® *HBV*-Monitor-FRT PCR kit should be transported at 2–8 °C for no longer than 5 days.

12. STABILITY AND STORAGE

All components of the PCR kit variant FRT, *HBV*-Q calibration kit and Negative Control (C-) are to be stored at the temperature from minus 24 to minus 16 °C when not in use. All components of the AmpliSens® *HBV*-Monitor-FRT PCR kit are stable until the expiry date stated on the label. The shelf life of the reagents before and after the first use is the same, unless otherwise stated.

NOTE: PCR-mix-1-FL *HBV* is to be kept away from light.

NOTE: Do not repeat freeze-thaw cycles more than twice for Positive Control-1-*HBV*, Positive Control-2-*HBV*, PIC1 *HBV*, PIC2 *HBV*, Internal Control STI-87. Store the abovementioned reagents at 2–8 °C for up to 6 months after thawing.

13. SPECIFICATIONS

13.1 Linear measurement range

The linear measurement range of AmpliSens® *HBV*-Monitor-FRT PCR kit is specified in the table below.

Table 6

Biological material	Extraction volume, µl	Nucleic acid extraction kit	Linear measurement range, IU/ml
Blood plasma	100	RIBO-prep NucliSENS easyMAG	150 – 100,000,000
	200	MAGNO-sorb	75 – 100,000,000
	1,000	MAGNO-sorb NucliSENS easyMAG	15 – 100,000,000

13.2 Analytical specificity

The analytical specificity of AmpliSens® *HBV*-Monitor-FRT PCR kit is ensured by the selection of specific primers and probes as well as stringent reaction conditions. The primers and probes were checked for possible homologies to all sequences published in the gene banks by sequence comparison analysis.

The analytical specificity is also ensured by the addition of the genomic RNA/DNA of the following organisms and viruses to the reaction: hepatitis A virus; hepatitis D virus; hepatitis C virus; human immunodeficiency virus; cytomegalovirus; Epstein-Barr virus; herpes simplex virus types 1 and 2; varicella-zoster virus; human herpes virus types 6 and 8; parvovirus B19; tick-borne encephalitis virus; West Nile encephalitis virus; adenovirus types 2, 3, and 7; *Escherichia coli*; *Staphylococcus aureus*; *Streptococcus pyogenes*; *Streptococcus agalactiae*; and *Homo sapiens*.

No cross-reactions were observed for the abovementioned organisms and viruses.

14. REFERENCES

- Handbook "Sampling, transportation, storage of clinical material for PCR diagnostics", developed by Federal State Institute of Science "Central Research Institute of Epidemiology" of Federal Service for Surveillance on Consumers' Rights Protection and Human Well-Being.
- Guidelines to AmpliSens® *HCV*-Monitor-FRT, AmpliSens® *HBV*-Monitor-FRT, and AmpliSens® *HDV*-Monitor-FRT PCR kits for quantitative detection of hepatitis C virus (*HCV*) RNA, hepatitis B virus (*HBV*) DNA and hepatitis D (*HDV*) RNA in the biological material by polymerase chain reaction (PCR) with real-time hybridization-fluorescence detection developed by Federal Budget Institute of Science "Central Research Institute for Epidemiology".

15. QUALITY CONTROL

In compliance with Federal Budget Institute of Science "Central Research Institute for Epidemiology" ISO 13485-Certified Quality Management System, each lot of AmpliSens® *HBV*-Monitor-FRT PCR kit has been tested against predetermined specifications to ensure consistent product quality.

List of Changes Made in the Instruction Manual		
VER	Location of changes	Essence of changes
03.12.10	Sampling and handling	Sentence «Collect blood samples into tubes with 3% EDTA solution (1 : 20) after overnight fasting» is changed into «Blood samples are taken after overnight fasting into the tube with EDTA solution as anticoagulant».
	Through the text	Corrections through the text MAGNO-sorb mention was deleted
	Stability and storage	Phrase about keeping away from light of PCR-mix-1-FL <i>HBV</i> is added
13.12.10	Content	DNA calibrators PIC1 <i>HBV</i> (C+1) and PIC2 <i>HBV</i> (C+2) are changed to DNA calibrators PIC1 <i>HBV</i> and PIC2 <i>HBV</i>
23.12.10 RT	Data analysis	Item was corrected in accordance with Russian instruction
	Stability and storage	<i>HBV</i> -Q calibration kit storage conditions are added
20.01.11 RT	Cover page	Phrase "For Professional Use Only" was added
	Intended use	The phrase "The results of PCR analysis are taken into account in complex diagnostics of disease" was added

VER	Location of changes	Essence of changes	
	Content	New sections "Working Conditions" and "Transportation" were added The "Explanation of Symbols" section was renamed to "Key to Symbols Used"	
	Stability and Storage	The information about the shelf life of open reagents was added	
	Key to Symbols Used	The explanation of symbols was corrected	
03.07.11 RT	Cover page, text	The name of Institute was changed to Federal Budget Institute of Science "Central Research Institute for Epidemiology"	
15.09.11 RT	8. PROTOCOL 8.1. DNA extraction	The information about using RIBO-prep kit was added	
02.12.11 VV	3.CONTENT	The mentioning about the PCR kit form that does not correspond to the specified catalogue number is deleted	
16.03.12 LA	13.1. Sensitivity	Column "Test material" (blood plasma) is added The name of DNA/RNA extraction kit is changed from RIBO sorb to RIBO-sorb-12	
20.06.12 LA	Cover page		
	16. Key to symbols used	Symbol [VD] was replaced by [RUQ] symbol	
	13.1. Sensitivity	Information related to MAGNO-sorb reagent kit was added to the table of analytical sensitivity	
	8.1. DNA extraction	Reference number of MAGNO-sorb reagent kit was added: K2-16-1000-CE Information about extraction with MAGNO-sorb is added	
	2. Principle of PCR detection	Paragraph revised in accordance to the Russian Instruction Manual and the template	
	4. Additional Requirements	DNase free filter tips were added to the list Corrections through the text	
23.06.14 ME	5 General Precautions	Section corrected with respect to the template	
	6. Sampling and Handling	Attention box was added about the blood samples storage Corrections through the text	
	8.1 DNA extraction	Section corrected, the use of the EM-plus kit was deleted, with respect to the template and Russian Instruction Manual section for NucliSENS easyMAG extraction was moved to the Guidelines, Calibration and calculation of coefficient B using <i>HBV-Q</i> calibration kit was added from Appendix 2.	
	8.2.1 Preparing tubes for PCR	Table 1 was added from Appendix 1 The names of positive controls of extraction were changed from "PCE 1" and "PCE 2" to "PCE-1" and "PCE-2"	
	9. Data analysis	Chapter rewritten	
	10. Troubleshooting	Section corrected with respect to the template and Russian Instruction Manual	
	13.1 Sensitivity	Table 5, the name of column 1 was changed from "Test material" to "Clinical material"	
	Throughout the text	Minor grammar and syntax corrections Tables were numbered	
	14.03.15 ME	Through the text	Misprints and inaccuracies was corrected
	08.05.15 ME	Text	Clinical material was changed to biological
1. Intended use		The phrase "The results of PCR analysis are taken into account in complex diagnostics of disease" was changed to "For research use only. Not for diagnostic procedures"	
3. Content		Compact Disk with software was added	
5. General precautions		Corrections in accordance with template	
9. Data analysis		The unit of measurement for coefficient B was deleted. In the table with results for controls criteria for the FAM channel were changed	
06.06.17 ME	Cover page	[REF] R-V5-P-M(RG,iQ,Mx,Dt)-CE was added	
	3. Content	Number of forms was changed from 1 to 2. It is specified that Form 1 is described in the Instruction Manual to AmpliSens® <i>HBV-Monitor-FRT</i> PCR kit [REF] TR-V5-P-M(RG,iQ,Mx,Dt)-CE. The content of form 2 was specified	
08.02.18 DV	Cover page	[REF] R-V5-MC(RG,iQ,Mx,Dt)-CE-B was added	
	3. Content	[REF] R-V5-MC(RG,iQ,Mx,Dt)-CE-B was added The third form was added – PCR kit variant FRT in bulk	
06.09.18 EM	3. Content	The colour of the reagent was specified	
30.06.20 KK	3. Content	The information about Form 1 [REF] TR-V5-P-M(RG,iQ,Mx,Dt)-CE was added	
17.07.20 MM	Throughout the text	The text formatting was changed	
	Footer	The phrase "For research use only. Not for diagnostic procedures" was added	
07.10.20 EM	Throughout the text	All the sections were updated according to the template	
	10. Troubleshooting	The information about the results not taken into account was deleted in case of absent the Ct values for Positive Controls and added in case of lower the concentration for IC	
16.11.20 MM	Content	The second form was deleted – PCR kit variant FRT in bulk	
	Footer	[REF] R-V5-MC(RG,iQ,Mx,Dt)-CE-B was deleted	
01.06.21 MM	Content	The form in bulk was added	
	Footer	The [REF] R-V5-MC(RG,iQ,Mx,Dt)-CE-B was added	

AmpliSens®



Federal Budget Institute of
Science "Central Research
Institute for Epidemiology"
3A Novogireevskaya Street
Moscow 111123 Russia