

Instruction Manual

KEY TO SYMBOLS USED

	Catalogue number		Contains sufficient for <n> tests
	Batch code		Use-by Date
	Research use only		Consult instructions for use
	Version		Keep away from sunlight
	Temperature limit	NCA	Negative control of amplification
	Manufacturer	C-	Negative control of extraction
	Date of manufacture	C+	Positive control of amplification
	Caution	IC	Internal control
		PCE	Positive Control of Extraction

1. INTENDED USE

AmpliSens® HCV / HBV / HIV-FRT PCR kit is an *in vitro* nucleic acid amplification test for simultaneous detection of hepatitis C virus RNA (HCV), hepatitis B virus DNA (HBV) and human immunodeficiency virus RNA (HIV) in the biological material using real-time hybridization-fluorescence detection of amplified products.

This kit can be used for analysis of both individual samples and several plasma samples combined in a minipool. For individual samples, the volume of plasma for nucleic acid extraction is 100, 200 or 1000 µl. For mini-pools, the recommended volume of sample is 1000 µl. The volume of each sample in mini-pool should be no less than 100 µl. The recommended number of plasma samples, combined to mini-pool, is 4–10 with the resulting volume of 1000 µl.

Extraction methods, volumes used for extraction, and the possibility of testing samples combined into mini-pools are specified in Table 1.

Table 1

Extraction method	Volumes of plasma used for RNA/DNA extraction, µl	Possibility of testing plasma mini-pools
NucliSENS easyMAG	100	no
	1,000	Up to 10 samples
MAGNO-sorb	200	no
	1,000	Up to 10 samples
QIASymphony Virus/Bacteria Midi Kit	1,000	Up to 10 samples
RIBO-sorb	100	no
RIBO-prep	100	no

NOTE: For research use only. Not for diagnostic procedures.

2. PRINCIPLE OF PCR DETECTION

Detection of hepatitis C virus RNA (HCV), hepatitis B virus DNA (HBV) and human immunodeficiency virus type 1 RNA (HIV-1) in biological material by using real-time hybridization-fluorescence detection includes the following stages: (1) RNA/DNA extraction from blood plasma together with Internal Control sample (IC) by means of QIASymphony SP and NucliSENS easyMAG automated systems or by means of manual extraction with RIBO-sorb, MAGNO-sorb or RIBO-prep kits; (2) reverse transcription of RNA and amplification of DNA/cDNA by using real-time hybridization-fluorescence detection.

The results of amplification are registered in the following fluorescence channels:

Table 2

Channel for fluorophore	FAM	JOE	ROX	Cy5	Cy5.5 ¹
DNA/cDNA-target	HCV cDNA	HIV-1 cDNA	HBV DNA	IC cDNA	HIV-2 cDNA
Target gene	5' UTR	Int/ LTR	gene C	Artificially synthesized sequence	LTR

¹ HIV-2 cDNA is not being detected when using variant FRT-4x.

3. CONTENT

AmpliSens® HCV / HBV / HIV-FRT PCR kit is produced in 4 forms:

variant FRT-A  R-V62-Q(RG,Dt)-CE adopted for the QIASymphony AS (QIAGEN, Germany) automatic module for reaction mixture preparation.

variant FRT  R-V62(RG,Dt)-CE;

variant FRT-4x  R-V50-4x(RG,iQ,Mx,Dt)-CE,

variant FRT in bulk²  R-V62(RG,Dt)-CE-B.

Variant FRT includes:

Reagent	Description	Volume, ml	Quantity
RT-G-mix-2	colorless clear liquid	0.015	4 tubes
RT-PCR-mix-1-FRT HCV / HBV / HIV1 / HIV2	clear liquid from colorless to light lilac colour	0.3	4 tubes
RT-PCR-mix-2-FL	colorless clear liquid	0.3	4 tubes
Polymerase (TaqF)	colorless clear liquid	0.03	4 tubes
TM-Revertase (MMIv)	colorless clear liquid	0.015	4 tubes
Buffer for elution	colorless clear liquid	1.2	4 tubes
Positive Control cDNA HCV / HBV / HIV (C+HCV / HBV / HIV)	colorless clear liquid	0.1	4 tubes
Negative Control (C-)*	colorless clear liquid	0.5	4 tubes
Positive Control HCV / HBV / HIV-rec**	colorless clear liquid	0.2	4 tubes
Internal Control STI-87-rec (IC)***	colorless clear liquid	0.5	4 tubes

* must be used in the extraction procedure as Negative Control of Extraction.

** must be used in the extraction procedure as Positive Control of Extraction.

*** add 10 µl of Internal Control STI-87-rec (IC) during the RNA/DNA extraction procedure directly to the sample/lysis mixture.

Variant FRT is intended for 100 reactions, including controls.

Variant FRT A includes:

Reagent	Description	Volume, ml	Quantity
RT-G-mix-2	colorless clear liquid	0.055	4 tubes
RT-PCR-mix-1-FRT HCV / HBV / HIV1 / HIV2	clear liquid from colorless to light lilac colour	0.35	4 tubes
RT-PCR-mix-2-FL	colorless clear liquid	0.35	4 tubes
Polymerase (TaqF)	colorless clear liquid	0.07	4 tubes
TM-Revertase (MMIv)	colorless clear liquid	0.055	4 tubes
Buffer for elution	colorless clear liquid	0.2	4 tubes
Positive Control cDNA HCV / HBV / HIV (C+HCV / HBV / HIV)	colorless clear liquid	0.1	4 tubes
Negative Control (C-)*	colorless clear liquid	1.2	4 tubes
Positive Control HCV / HBV / HIV-rec**	colorless clear liquid	1.2	4 tubes
Internal Control STI-87-rec (IC)***	colorless clear liquid	0.5	4 tubes

* must be used in the extraction procedure as Negative Control of Extraction.

** must be used in the extraction procedure as Positive Control of Extraction.

*** add 10 µl of Internal Control STI-87-rec (IC) during the RNA/DNA extraction procedure

Variant FRT A is intended for 100 reactions, including controls.

Variant FRT-4x includes:

Reagent	Description	Volume, ml	Quantity
RT-G-mix-2	colorless clear liquid	0.015	4 tubes
RT-PCR-mix-1-FRT HCV / HBV / HIV	clear liquid from colorless to light lilac colour	0.3	4 tubes
RT-PCR-mix-2-FL	colorless clear liquid	0.3	4 tubes
Polymerase (TaqF)	colorless clear liquid	0.03	4 tubes
TM-Revertase (MMIv)	colorless clear liquid	0.015	4 tubes
Buffer for elution	colorless clear liquid	1.2	4 tubes
Positive Control cDNA HCV / HBV / HIV (C+HCV / HBV / HIV)	colorless clear liquid	0.1	4 tubes
Negative Control (C-)*	colorless clear liquid	0.5	4 tubes
Positive Control HCV / HBV / HIV-rec**	colorless clear liquid	0.2	4 tubes
Internal Control STI-87-rec (IC)***	colorless clear liquid	0.5	4 tubes

* must be used in the extraction procedure as Negative Control of Extraction.

** must be used in the extraction procedure as Positive Control of Extraction.

*** add 10 µl of Internal Control STI-87-rec (IC) during the RNA/DNA extraction procedure directly to the sample/lysis mixture

Variant FRT-4x is intended for 100 reactions, including controls.

² In bulk form contains unlabeled tubes. Tubes with identical reagent are packed in one bag with label.

4. ADDITIONAL REQUIREMENTS

- DNA/RNA extraction kit or RNA/DNA extraction automatic station.
- Disposable powder-free gloves and laboratory coat.
- Pipettes (adjustable).
- Sterile RNase-free pipette tips with aerosol filters (up to 200 µl).
- Tube racks.
- Vortex mixer.
- Desktop centrifuge with rotor for 2-ml reaction tubes.
- PCR box.
- For variant FRT-A: QIASymphony SP + AS (QIAGEN, Germany) automatic station for RNA/DNA extraction and reaction mixture preparation and consumable materials for QIASymphony AS (QIAGEN, Germany) automatic module (filter tips up to 1,500, 200 and 50 µl).
- For variant FRT-A and variant FRT: Real-time instruments (for example, Rotor-Gene 3000 or Rotor-Gene 6000 (Corbett Research, Australia); CFX96 (Bio-Rad, USA) or equivalent) with possibility of simultaneous detection in five channels (corresponding to FAM, JOE, ROX, Cy5, Cy 5.5 fluorescent dyes).
- For variant FRT-4x: Real-time instruments (for example, Rotor-Gene 3000 or Rotor-Gene 6000 (Corbett Research, Australia); CFX96 (Bio-Rad, USA) or equivalent) with possibility of simultaneous detection in four or more channels (corresponding to FAM, JOE, ROX, Cy5 fluorescent dyes).
- Disposable polypropylene tubes for PCR (0.1- or 0.2-ml)
 - 0.2-ml PCR tubes or strips with domed caps (for example, Axygen, USA) if a plate-type instrument is used;
 - 0.2-ml PCR tubes (flat caps, nonstriped) (for example, Axygen, USA) for 36-well rotor or 0.1-ml tubes (Corbett Research, Australia) for 72-well rotor if a rotor-type instrument is used.
- Refrigerator at the temperature from 2 to 8 °C.
- Deep-freezer at the temperature from minus 24 to minus 16 °C.
- Reservoir for used tips.

5. GENERAL PRECAUTIONS

The user should always pay attention to the following:

- Use sterile pipette tips with aerosol filters and use a new tip for every procedure.
- Store all extracted positive material (specimens, controls and amplicons) away from all other reagents and add it to the reaction mix in a distantly separated facility.
- Thaw all components thoroughly at room temperature before starting an assay.
- When thawed, mix the components and centrifuge briefly.
- Use disposable protective gloves and laboratory cloths, and protect eyes while samples and reagents handling. Thoroughly wash hands afterwards.
- Do not eat, drink, smoke, apply cosmetics, or handle contact lenses in laboratory work areas.
- Do not use a kit after its expiration date.
- Dispose of all specimens and unused reagents in accordance with local regulations.
- Samples should be considered potentially infectious and handled in biological cabinet in compliance with appropriate biosafety practices.
- Clean and disinfect all samples or reagents spills using a disinfectant, such as 0.5 % sodium hypochlorite or another suitable disinfectant.
- Avoid inhalation of vapors, samples and reagents contact with the skin, eyes, and mucous membranes. Harmful if swallowed. If these solutions come into contact, rinse the injured area immediately with water and seek medical advice if necessary.
- Safety Data Sheets (SDS) are available on request.
- Use of this product should be limited to personnel trained in DNA amplification techniques.
- Workflow in the laboratory must be one-directional, beginning in the Extraction Area and moving to the Amplification and Detection Area. Do not return samples, equipment and reagents in the area where the previous step was performed.



Some components of this kit contain sodium azide as a preservative. Do not use metal tubing for reagent transfer.

6. SAMPLING AND HANDLING

NOTE: Obtaining samples of biological materials for PCR-analysis, transportation and storage are described in the manufacturer's handbook [1]. It is recommended that this handbook is read before starting work.

Analysis is performed with peripheral blood plasma. Blood is obtained after overnight fasting to tube with EDTA solution as anticoagulant. Turn the closed tube several times. Plasma should be transferred to a new tube within 6 h after bleeding. This is performed by centrifuging the tube with blood at 800–1600 g for 20 min. Plasma can be stored for no more than 5 days at 2–8 °C and for a long time at the temperature not more than minus 68 °C.

In some cases, blood serum can be used. The analytical sensitivity of the reagent kit for this material is the same, but the clinical sensitivity may greatly decrease because of sedimentation of viral particles as a result of retraction of the clot. Store blood serum for no more than 5 days at 2–8 °C and for a long time at the temperature not more than minus 68 °C.

This kit can be used for analyzing both individual samples and several plasma samples combined into a mini-pool. The recommended number of plasma samples for pooling is 4–10 (the resulting volume, 1000 µl).

7. WORKING CONDITIONS

AmpliSens® HCV / HBV / HIV-FRT PCR kit should be used at 18–25 °C.

8. PROTOCOL

8.1. RNA extraction

It's recommended to use the following nucleic acid extraction kits:

- RIBO-sorb,
- RIBO-prep,
- MAGNO-sorb.

QIASymphony SP and NucliSENS easyMAG automated system can be used.

If using RIBO-prep kit, extract the RNA/DNA according to the manufacturer's protocol taking into account next additions and improvements:

- It is allowed to mix the **Solution for lysis** and **Internal Control STI-87-rec (IC)** in a separate sterile vial (300 µl of **Solution for lysis** and 10 µl of **Internal Control STI-87-rec (IC)** per sample) and then transfer 300 µl of mixture to each prepared 1.5-ml tube to simplify the extraction procedure in case of great quantity of samples.
- To prepare the Positive Control of extraction, **PCE**, add 100 µl of the **Positive Control HCV / HBV / HIV-rec** sample to a tube containing **Solution for Lysis**.

NOTE:

If using RIBO-sorb kit, extract the RNA/DNA according to the manufacturer's protocol taking into account next additions and improvements:

- Add 10 µl of **Internal Control STI-87-rec (IC)** to each tube and then add 450 µl of **Lysis Solution**.
- It is allowed to mix the **Lysis Solution** and **Internal Control STI-87-rec (IC)** in a separate sterile vial (450 µl of **Lysis Solution** and 10 µl of **Internal Control STI-87-rec (IC)** per sample) and then transfer 450 µl of mixture to each prepared 1.5-ml tube to simplify the extraction procedure in case of great quantity of samples.
- To prepare the Positive Control of extraction, **PCE**, add 100 µl of the **Positive Control HCV / HBV / HIV-rec** sample to a tube containing **Lysis Solution**.

NOTE:

If using the MAGNO-sorb kit extract the RNA/DNA according to the manufacturer's protocol taking into account next additions and improvements:

- In case of RNA/DNA extraction from blood plasma sample of 1000 µl, the volume of the **Internal Control STI-87-rec (IC)** required for 24-tube panel is 0.28 ml. In case of other panels and RNA/DNA extraction from blood plasma sample of 200 µl see MAGNO-sorb instruction manual.

NOTE:

- To prepare the Positive Control of extraction, **PCE**, add 100 µl of the **Positive Control HCV / HBV / HIV-rec** sample to a tube containing **Lysis Solution MAGNO-sorb**.
- To prepare the Negative Control of extraction, **C-**, add 100 µl of the **Negative Control (C-)** sample to a tube containing **Lysis Solution MAGNO-sorb**.
- The volume of **Buffer for elution** required for extraction from both 1000 and 200 µl of blood plasma samples is 70 µl.

NOTE: **If using the QIASymphony SP automated system** extract the RNA/DNA according to the Guidelines [2].

8.2. Reverse transcription

Total reaction volume is 50 µl, the volume of RNA/DNA sample is 30 µl.

NOTE: RNase-free and Dnase-free plastic ware should be used only.

8.2.1 Preparing tubes for RT-PCR

Variant FRT

Tubes type depends on the thermocycler used. Use filter tips for transferring the reagents and samples into the tubes.

NOTE: All components of the reaction mixture should be mixed immediately before use.

1. Before starting work, thaw and thoroughly vortex all reagents of the kit. Sediment the drops from the tubes' caps.
2. Take the required number of PCR-tubes for clinical and control samples (2 controls of the extraction stage and 2 controls of amplification stage). The type of tubes, strips or plates depends on the real-time PCR instrument used for analysis.
3. Prepare the reaction mixture (see Table 3) for the required number of reactions including the analysis of test and control samples.

Table 3

Reaction mixture preparation						
		Reagent volume per specified number of PCR reactions with extra volume				
Reagent volume per 1 reaction, µl		10.00	10.00	1.00	0.50	0.50
Number of clinical samples	Number of PCR reactions ³	RT-PCR-mix-1-FRT HCV / HBV / HIV1 / HIV2	RT-PCR-mix-2-FL	Polymerase (TaqF)	RT-G-mix-2	TM-Revertase (MMIv)
6	10	110	110	11	5	5
8	12	130	130	13	6	6
10	14	150	150	15	7	7
12	16	170	170	17	8	8
14	18	190	190	19	9	9
16	20	210	210	21	10	10
18	22	230	230	23	11	11
20	24	250	250	25	12	12
22 ⁴	26	Whole tube content	Whole tube content	Whole tube content	15	Whole tube content

4. For analysis of 22 clinical samples (24 samples after extraction procedure) the following method is recommended: transfer the whole volume of **RT-PCR-mix-2-FL** (0.3 ml), the whole volume of **polymerase (TaqF)** (0.03 ml), the whole volume of **TM-Revertase (MMIv)** (0.015 ml) and the whole volume of **RT-G-mix-2** (0.015 ml) into the tube with **RT-PCR-mix-1-FRT HCV / HBV / HIV1 / HIV2** (0.3 ml). Thoroughly mix by vortexing and sediment the drops from the tubes caps.

5. Add 20 µl of prepared reaction mixture into each PCR-tube. Discard the remaining mixture.
6. Using filter tips add 30 µl of RNA/DNA-sample to reaction mixture. Thoroughly mix the content by pipetting avoiding foaming.

³ Number of clinical samples + 2 controls for RNA extraction stage + 2 controls of RT-PCR (N+4, N – number of clinical samples).

⁴ Typical extraction from 24 samples using NucliSENS easyMAG automated system or MAGNO-sorb nucleic acids extraction kit.

7. Carry out the control amplification reactions:	
C+	— Add 30 µl of Positive Control cDNA HCV / HBV / HIV (C+_{HCV/ HBV/ HIV}) to the tube labeled C+.
NCA	— Add 30 µl of Buffer for elution to the tube labeled NCA.
C-	— Add 30 µl of the sample extracted from the Negative Control (C-) reagent to the tube labeled C- (Negative control of Extraction).
PCE	— Add 30 µl of the sample extracted from the Positive Control HCV / HBV / HIV-rec reagent to the tube labeled PCE (Positive control of Extraction).

NOTE: Tubes with control samples are to be mixed by pipetting avoiding foaming.

Variant FRT-A

Variant FRT-A is adopted for the QIASymphony AS (QIAGEN, Germany) automatic module for reaction mixture preparation. Preparing the reaction mixtures, adding the RNA/DNA samples and controls into the tubes with reaction mixture are performed automatically. The AmpliSens[®] RG protocol, preset in the instrument, should be used for **AmpliSens[®] HCV / HBV / HIV-FRT** PCR kit. The work with the QIASymphony AS automatic module is to be carried out in accordance with the instrument manual.

NOTE: Thaw and thoroughly mix on vortex all reagents of the kit before work. Sediment the drops from the tubes' caps. Take the required number of tubes for amplification of clinical and control samples.

NOTE: The control amplification reactions (C+ and NCA) are allowed to be omitted in case of reaction mixtures preparation with the use of automatic module.

Variant FRT-4x

Tubes type depends on the thermocycler used. Use filter tips for transferring the reagents and samples into the tubes.

NOTE: All components of the reaction mixture should be mixed immediately before use.

- Before starting work, thaw and thoroughly vortex all reagents of the kit. Sediment the drops from the tubes' caps.
- Take the required number of PCR-tubes for clinical and control samples (2 controls of the extraction stage and 2 controls of amplification stage). The type of tubes, strips and plates depends on the real-time PCR instrument used for analysis.
- Prepare the reaction mixture (see Table 4) for the required number of reactions including the analysis of test and control samples.

Table 4

Reaction mixture preparation						
Reagent volume per specified number of PCR reactions with extra volume						
Reagent volume for 1 reaction, µl		10.00	10.00	1.00	0.50	0.50
Number of clinical samples	Number of PCR reactions ⁵	RT-PCR-mix-1-FRT HCV / HBV / HIV	RT-PCR-mix-2-FL	Polymerase (TaqF)	RT-G-mix-2	TM-Revertase (MMIV)
6	10	110	110	11	5	5
8	12	130	130	13	6	6
10	14	150	150	15	7	7
12	16	170	170	17	8	8
14	18	190	190	19	9	9
16	20	210	210	21	10	10
18	22	230	230	23	11	11
20	24	250	250	25	12	12
22 ⁶	26	Whole tube content	Whole tube content	Whole tube content	15	Whole tube content

4. For analysis of 22 clinical samples (24 samples after extraction procedure) the following method is recommended: transfer the whole volume of **RT-PCR-mix-2-FL** (0.3 ml), the whole volume of **polymerase (TaqF)** (0.03 ml), the whole volume of **TM-Revertase (MMIV)** (0.015 ml) and the whole volume of **RT-G-mix-2** (0.015 ml) into the tube with **RT-PCR-mix-1-FRT HCV / HBV / HIV** (0.3 ml). Thoroughly mix by vortexing and sediment the drops from the tubes caps.

5. Add 20 µl of prepared reaction mixture into each PCR-tube. Discard the remaining mixture.

6. Using filter tips add 30 µl of RNA/DNA-sample to reaction mixture. Mix the content by pipetting avoiding foaming.

7. Carry out **control amplification reactions**:

C+ — Add 30 µl of **Positive Control cDNA HCV / HBV / HIV (C+_{HCV/ HBV/ HIV})** to the tube labeled C+.

NCA — Add 30 µl of **Buffer for elution** to the tube labeled NCA.

C- — Add 30 µl of the **sample extracted from the Negative Control (C-)** reagent to the tube labeled C- (Negative control of Extraction).

PCE — Add 30 µl of the **sample extracted from the Positive Control HCV / HBV / HIV-rec** reagent to the tube labeled PCE (Positive control of Extraction).

NOTE: Tubes with control samples are to be mixed by pipetting avoiding foaming.

8.2.2. Amplification

- Create a temperature profile on your instrument as follows (see Tables 5, 6):

Table 5

AmpliSens [®] HCV / HBV / HIV amplification program for rotor-type instruments ⁷				
Step	Temperature, °C	Time	Fluorescence detection	Cycles
1	50	20 min	—	1
2	95	15 min	—	1
3	95	20 s	—	4
	46	40 s	—	
4	95	5 s	—	42
	60	40 s	—	
	45	30 s	FAM, JOE, ROX, Cy5, Cy5.5	

NOTE: When using **variant FRT-4x** the channel for the Cy5.5 fluorophore is not used.

⁵ Number of clinical samples + 2 controls for RNA extraction stage + 2 controls of RT-PCR (N+4, N – number of clinical samples).

⁶ Typical extraction from 24 samples using NucliSENS easyMAG automated system or MAGNO-sorb nucleic acids extraction kit.

⁷ For example, Rotor-Gene 3000 and Rotor-Gene 6000 (Corbett Research, Australia).

Table 6
AmpliSens[®] HCV / HBV / HIV amplification program for plate-type instruments⁸

Step	Temperature, °C	Time	Fluorescence detection	Cycles
1	50	20 min	—	1
2	95	15 min	—	1
3	95	20 s	—	4
	46	40 s	—	
4	95	5 s	—	42
	60	40 s	—	
	40	40 s	FAM, JOE, ROX, Cy5, Cy5.5	

NOTE: When using **variant FRT-4x** the channel for the Cy5.5 fluorophore is not used.

2. Adjust the fluorescence channel sensitivity according to the *Important Product Information Bulletin* and Guidelines [2].

3. Insert tubes into the reaction module of the device.

It is recommended to sediment drops from walls of tubes by short centrifugation

NOTE: (1–3 s) before placing them into the plate-type instrument.

4. Run the amplification program with fluorescence detection.

5. Analyze results after the amplification program is completed.

9. DATA ANALYSIS

The results are interpreted by the real-time PCR instrument software:

- HCV cDNA amplification product is detected in the channel for the FAM fluorophore.
- HIV-1 cDNA amplification product is detected in the channel for the JOE fluorophore.
- HBV DNA amplification product is detected in the channel for the ROX fluorophore.
- IC cDNA amplification product is detected in the channel for the Cy5 fluorophore.
- HIV-2 cDNA amplification product is detected in the channel for the Cy5.5 fluorophore.

The results are interpreted by the crossing (or not crossing) the threshold line by the fluorescence curve that indicates the presence (or absence) of the threshold Ct value in corresponding column in the results table.

The result of amplification is considered *positive* if the fluorescence S-shaped curve (typical for real-time PCR) crosses once the threshold line in the fluorescence reliable growth area; and the threshold value (Ct or Cp) for respective channel is less than the Ct value specified in the *Important Product Information Bulletin*.

The result of amplification is considered *negative* if the fluorescence curve is not S-shaped or does not cross the threshold line (Ct or Cp values are absent).

Otherwise, the result of amplification is considered *equivocal*.

For details, see Instruction Manuals for appropriate instruments and Guidelines [2].

The sample is considered **positive** for **HBV DNA, HCV RNA, HIV-1 RNA and HIV-2 RNA** if the Ct value defined in respective channel is less than the boundary Ct value specified in the *Important Product Information Bulletin* and the Ct value in the Cy5 channel is less than the specified boundary Ct value.

The sample is considered **negative** for **HBV DNA, HCV RNA, HIV-1 RNA and HIV-2 RNA** if the Ct value defined in respective channel is absent or greater than the boundary Ct value specified in the *Important Product Information Bulletin* and Ct value in Cy5 channel is less than the specified boundary Ct value.

NOTE: Boundary Ct values are specified in the *Important Product Information Bulletin* enclosed in the PCR kit.

NOTE: When using **variant FRT-4x** the channel for the Cy 5.5 fluorophore is not used and the analysis for this channel is not performed.

The results of the analysis are considered **reliable only if the results for control samples are correct** (see Table 7).

Table 7

Results for controls					
Control	Ct value in the channel for the fluorophore				
	FAM	JOE	ROX	Cy5	Cy5.5
C-	absent	absent	absent	≤ boundary value (positive)	absent
PCE	≤ boundary value (positive)	≤ boundary value (positive)	≤ boundary value (positive)	≤ boundary value (positive)	≤ boundary value (positive)
NCA	absent	absent	absent	absent	absent
C+	≤ boundary value (positive)	≤ boundary value (positive)	≤ boundary value (positive)	≤ boundary value (positive)	≤ boundary value (positive)

NOTE: Boundary Ct values are specified in the *Important Product Information Bulletin* enclosed in the PCR kit.

10. TROUBLESHOOTING

- If the Ct value in the channels for FAM, JOE, ROX or Cy5.5 fluorophores is greater than the boundary Ct value and the Ct value in the Cy5 channel is also greater than the boundary Ct value, PCR should be repeated for this sample starting from the RNA/DNA extraction stage.
- If the Ct value for the Positive Control of Extraction (PCE) or Positive Control of Amplification (C+) in any channel is absent or greater than the specified boundary Ct value, PCR should be repeated for all the samples starting from the RNA/DNA extraction stage.
- If the Ct value is present for the Negative Control of Extraction (C-) in the channels FAM, JOE, ROX, Cy5.5 fluorophores and/or Negative Control of Amplification (NCA) in any channel, PCR should be repeated for all positive samples starting from the RNA/DNA extraction stage.

11. TRANSPORTATION

AmpliSens[®] HCV / HBV / HIV-FRT PCR kit should be transported at 2–8 °C for no longer than 5 days.

⁸ For example, CFX96 (Bio-Rad, USA).

12. STABILITY AND STORAGE

All components of the **AmpliSens® HCV / HBV / HIV-FRT** PCR kit are to be stored at temperature from minus 24 to minus 16 °C when not in use. All components of the **AmpliSens® HCV / HBV / HIV-FRT** PCR kit are to be stable until labeled expiration date. The shelf life of reagents before and after the first use is the same, unless otherwise stated.

NOTE: RT-PCR-mix-1-FRT **HCV / HBV / HIV1 / HIV2** and RT-PCR-mix-1-FRT **HCV / HBV / HIV** are to be kept away from light.

13. SPECIFICATIONS

13.1. Sensitivity

Analytical sensitivity of **AmpliSens® HCV / HBV / HIV-FRT** PCR kit is presented in the table below.

Volume of sample for extraction, µl	Method of extraction	Analytical sensitivity			
		HCV, IU/ml	HBV, IU/ml	HIV-1, copies/ml	HIV-2*, copies/ml
100	RIBO-sorb, RIBO-prep, NucliSENS easyMAG	100	50	200	600
200	MAGNO-sorb	50	25	100	300
1000	MAGNO-sorb, NucliSENS easyMAG, QIASymphony Virus/Bacteria Midi Kit	10	5	20	60

NOTE: * When using variant FRT-4x the HIV-2 RNA is not being detected.

The claimed analytical features of **AmpliSens® HCV / HBV / HIV-FRT** PCR kit are guaranteed only when additional reagents kits RIBO-sorb, RIBO-prep, MAGNO-sorb (manufactured by Federal Budget Institute of Science "Central Research Institute for Epidemiology"), or QIASymphony SP or NucliSENS easyMAG automated system kit are used.

NOTE:

13.2. Specificity

The analytical specificity of **AmpliSens® HCV / HBV / HIV-FRT** PCR kit is ensured by selection of specific primers and probes and stringent reaction conditions. The primers and probes were tested for possible homologues to all sequences deposited in gene banks by sequence comparison analysis as well as with genomic DNA/RNA of the following organisms and viruses: hepatitis A virus; hepatitis D virus; cytomegalovirus; Epstein-Barr virus; herpes simplex virus types 1 and 2; *Varicella-Zoster virus*; human herpes virus types 6 and 8; parvovirus B19; tick-borne encephalitis virus; *West Nile virus*; adenovirus types 2, 3, and 7; *Escherichia coli*, *Staphylococcus aureus*; *Streptococcus pyogenes*, *S.agalactiae*; and *Homo sapiens*.

14. REFERENCES

- Handbook "Sampling, Transportation, and Storage of Clinical Material for PCR diagnostics", developed by Federal Budget Institute of Science "Central Research Institute for Epidemiology" of Federal Service for Surveillance on Consumers' Rights Protection and Human Well-Being.
- Guidelines to **AmpliSens® HCV / HBV / HIV-FRT** PCR kit for simultaneous detection of hepatitis C virus RNA (*HCV*), hepatitis B virus DNA (*HBV*) and human immunodeficiency virus RNA (*HIV*) in the biological material by the polymerase chain reaction (PCR) with real-time hybridization-fluorescence detection developed by Federal Budget Institute of Science "Central Research Institute for Epidemiology".

15. QUALITY CONTROL

In compliance with Federal Budget Institute of Science "Central Research Institute for Epidemiology" ISO 13485-Certified Quality Management System, each lot of **AmpliSens® HCV / HBV / HIV-FRT** PCR kit has been tested against predetermined specifications to ensure consistent product quality.

List of Changes Made in the Instruction Manual

VER	Location of changes	Essence of changes
25.08.10	Page footer	Reference number is changed from R-V50-4x(RG,iQ,Mx,Dt) to R-V50-4x(RG,iQ,Mx,Dt)-CE
	3. Content	
08.07.11 LA	Cover page	The phrase "For Professional Use Only" was added
	Intended use	The phrase "The results of PCR analysis are taken into account in complex diagnostics of disease" was added
	Content	New sections "Working Conditions" and "Transportation" were added The "Explanation of Symbols" section was renamed to "Key to Symbols Used"
	Stability and Storage	The information about the shelf life of reagents before and after the first use was added Information that PCR-mix-1-FRT is to be kept away from light was added
	Key to Symbols Used	The explanation of symbols was corrected
	Cover page, text	The name of Institute was changed to Federal Budget Institute of Science "Central Research Institute for Epidemiology"
10.09.11 RT	Text	Corrections through the text Mention about MAGNO-sorb kit was deleted QIASymphony automated system was added Channel for detection of <i>HBV</i> DNA was changed from ROX/TexasRed/Orange to ROX/Orange PCR kit variant FRT was changed into PCR kit variant FRT-4x
	Content	The volumes of Negative Control (C-) and Positive Control <i>HCV / HBV / HIV-rec</i> are changed from 1.2 ml and 0.3 ml into 0.5 ml and 0.2 ml respectively
	Sampling and handling	Procedure of preparing plasma samples was corrected Information about blood serum was added
	8. Protocol	
	8.2. Preparing the reverse transcription and PCR	The volume of added Buffer for elution was changed from 10 µl into 30 µl
	8.3. Reverse transcription and	Table with amplification program for plate-type instruments was deleted

	amplification	
	9. Data analysis	Item was corrected (information about interpretation of results was added) Table "Results for controls" was changed
	10. Troubleshooting	Item was changed
	13. Specifications	Information about genomic DNA/RNA of different organisms and viruses was added
	13.2. Specificity	
10.07.12 IVI	Title page, Key to symbols used	Symbol IVD <i>in vitro</i> diagnostic medical device was changed to RUO research use only
17.07.13 ME	Content	The form in bulk was added
	References	The link for guidelines was added
	Page footer	REF R-V50-4x(RG,iQ,Mx,Dt)-CE-B was added
04.03.14 ME	8.1. RNA/DNA extraction	MAGNO-sorb kit was added. Additions and improvements of using RIBO-sorb and RIBO-prep kits were added
	8.2. Preparing reverse transcription and PCR	Table 1 was added from Appendix. The numeration of the tables was corrected Table 3 "Program AmpliSens <i>HBV / HCV / HIV</i> for plate-type instruments" was added
	14. References	The references were corrected
10.02.15 ME	Footer	REF R-V62(RG,Dt)-CE was added
	Content	PCR kit variant FRT was added
	Principle of PCR detection, Additional requirements, Protocol, Data analysis, Troubleshooting, Stability and storage, Specifications	Additions about using the PCR kit variant FRT was added
15.05.15 PM	Text	Clinical material was changed to biological
	1. Intended use	The phrase "The results of PCR analysis are taken into account in complex diagnostics of disease" was changed to "For research use only. Not for diagnostic procedures"
	13.2. Specificity	The phrase "Specificity of AmpliSens® <i>HCV / HBV / HIV-FRT</i> PCR kit was confirmed in laboratory clinical trials" was deleted
15.07.15 ME	Footer	REF R-V62-Q(RG,Dt)-CE was added
	Content	PCR kit variant FRT-A was added
	5. Additional requirements, 8. Protocol	Additions about using the PCR kit variant FRT-A was added
24.04.18 PM	3. Content	The colour of reagents was specified
30.07.20 KK	Through the text	The text formatting was changed
	Footer	The phrase "For research use only. Not for diagnostic procedures" was added
02.11.20 MM	2. Principle of PCR detection	The table with targets was added
16.11.20 KK	Footer	The REF R-V50-4x(RG,iQ,Mx,Dt)-CE-B was deleted
	3. Content	The form in bulk was deleted
	Footer	The REF R-V62(RG,Dt)-CE-B was added
04.08.23 BA	Through the text	The reference numbers of nucleic acid extraction kits were deleted
	3. Content	Variant FRT in bulk was added

AmpliSens®



Federal Budget Institute of Science "Central Research Institute for Epidemiology"
3A Novogireevskaya Street
Moscow 111123 Russia