



RUO

For Professional Use Only

Application note

Detection of Chlamydia trachomatis using
**eSens Ureaplasma/M.genitalium/M.hominis
QL PCR kit, REF ES3048A**

for Research Use Only. Not for use in diagnostic procedures.

PURPOSE OF THIS APPLICATION NOTE

The purpose of this application note is to provide an **extended guidance** for utilizing the eSens Ureaplasma/M.genitalium/M.hominis QL PCR kit. This document aims to assist laboratory professionals and researchers in accurately detecting the DNA of **not only** Ureaplasma spp. (U. parvum and U. urealyticum), Mycoplasma genitalium, and Mycoplasma hominis in various clinical samples, but **it also introduces the supplementary feature of the kit that enables the detection of Chlamydia trachomatis DNA for research use only (RUO)**, despite it not being CE IVD validated for this pathogen.

This application note provides specific supplementary data **highlighted by red color** to facilitate the use of this feature.

Application note

Detection of Chlamydia trachomatis using eSens Ureaplasma/M.genitalium/M.hominis QL PCR kit, REF ES3048A for RUO; based on IFU Version 01_04/2022

Ecoli Dx, s.r.o., Purkyňova 74/2, 110 00 Praha 1, The Czech Republic

PRINCIPLE OF PCR DETECTION

Table 1

Channel for fluorophore	FAM	JOE	ROX	Cy5	Cy5.5
DNA-target	<i>Chlamydia trachomatis (RUO)</i>	<i>Ureaplasma sp.</i>	<i>Mycoplasma genitalium</i>	Internal Control-FL	<i>Mycoplasma hominis</i>
Target gene	<i>cryptic plasmid</i>	<i>UreC</i>	<i>gyrB</i> gene	genetically engineered construction	16s rRNA gene

PROTOCOL

Amplification

Fluorescent signal is detected in the channels for the FAM, JOE, ROX, Cy5 and Cy5.5 fluorophores.

Instrument Settings

Test settings for rotor-type instruments

Channel	Calibrate/Gain Optimisation	Threshold	Dynamic tube	Slope Correct	More Settings/ Outlier Removal
FAM/Green	from 5 FI to 10 FI	0.1	On	Off	0%
JOE/Yellow	from 4 FI to 8 FI	0.1	On	Off	5%
ROX/Orange	from 4 FI to 8 FI	0.1	On	Off	5%
Cy5/Red	from 4 FI to 8 FI	0.07	On	On	5%
Cy5.5/Crimson	from 4 FI to 8 FI	0.1	On	On	20%

Test settings for plate-type instruments

Channel	Threshold
FAM, HEX, ROX, Cy5, Cy5.5	For each channel in Log Scale set the threshold line at the level of 10-20 % of maximum fluorescence obtained for the Positive Control of Amplification (C+) in the last amplification cycle.

DATA ANALYSIS

Analysis of results is performed by the software of the real-time PCR instrument used by measuring fluorescence signal accumulation in five channels:

Application note

Detection of *Chlamydia trachomatis* using eSens *Ureaplasma/M.genitalium/M.hominis* QL PCR kit, REF ES3048A for RUO; based on IFU Version 01_04/2022

Ecoli Dx, s.r.o., Purkyňova 74/2, 110 00 Praha 1, The Czech Republic

- The signal of the **Chlamydia trachomatis DNA** amplification product is detected in the channel for the **FAM** fluorophore

Results are interpreted by the crossing (or not-crossing) the fluorescence curve with the threshold line set at the specific level that corresponds to the presence (or absence) of a Ct value of the DNA sample in the corresponding column of the results grid.

Principle of interpretation is the following:

- **Chlamydia trachomatis DNA is detected** if the Ct value is determined in the results grid in the channel for the **FAM** fluorophore. Moreover, the fluorescence curve of the sample should cross the threshold line in the area of typical exponential growth of fluorescence.
- **Chlamydia trachomatis, Ureaplasma sp. (U.parvum and U.urealyticum), Mycoplasma genitalium, and Mycoplasma hominis DNA are not detected** in a sample if the Ct value is not determined (absent) (fluorescence curve does not cross the threshold line) in the channels for **FAM**, **JOE**, **ROX** and **Cy5.5** fluorophores whereas the Ct value determined in the channel for the **Cy5** fluorophore is less than the boundary Ct value.
- The result is **invalid** if the Ct value is not determined (absent) in the channel for **Cy5** fluorophore or greater than the specified boundary Ct value, whereas the Ct value in the channel for the **FAM**, **JOE**, **ROX** and **Cy5.5** fluorophores is not determined (absent) or greater than the specified boundary Ct value. In such cases, the PCR analysis should be repeated starting from the DNA extraction stage

Table 3

Results for controls

Control	Stage for control	Ct value in the channel for fluorophore	
		FAM, JOE, ROX and Cy5.5	Cy5
C-	DNA extraction	Absent	<boundary value
NCA	PCR	Absent	Absent
C+	PCR	<boundary value	<boundary value

Table 4

Boundary Ct values

Sample	Rotor-type instrument					Plate-type instrument				
	Channel for fluorophore									
	FAM	JOE	ROX	Cy5.5	Cy5	FAM	JOE	ROX	Cy5.5	Cy5
C+	35	35	35	35	33	38	38	38	38	36
C-	Ct is absent				33	Ct is absent				36
NCA	Ct is absent					Ct is absent				
Test samples	-	-	-	-	33	-	-	-	-	36

Application note

Detection of Chlamydia trachomatis using eSens Ureaplasma/M.genitalium/M.hominis QL PCR kit, REF ES3048A for RUO; based on IFU Version 01_04/2022

Ecoli Dx, s.r.o., Purkyňova 74/2, 110 00 Praha 1, The Czech Republic

SPECIFICATIONS

Sensitivity

The analytical sensitivity for *Chlamydia trachomatis*, *Ureaplasma* spp., *Mycoplasma genitalium*, and *Mycoplasma hominis* is not less than 5×10^2 genome equivalents per 1 ml of sample (GE/ml).

Application note

Detection of *Chlamydia trachomatis* using eSens *Ureaplasma/M.genitalium/M.hominis* QL PCR kit, REF ES3048A for RUO; based on IFU Version 01_04/2022

Ecoli Dx, s.r.o., Purkyňova 74/2, 110 00 Praha 1, The Czech Republic