

# Pipette Filter Tips

- Sterile and clean of contaminants (EO sterilization)
- Filter preventing contamination of your sample and the pipette cone
- Suitable for DNA/RNA extraction and PCR setup
- Compatible with Finn pipettes and other brands
- Convenient lid opening using one hand  
– 2 positions of opening: half (vertical), full (upside down)



Catalog number	Product name	Finntip alternative	Pipette cone size & color-coding
Actual max. volume	Label on box	Finnpipette compatibility	Tip length
			
BK-PT-10SF	Sterile tip with filter 10µl, 10x96 PCS/box	Finntip Filter 10	micro
10µl	10µl FILTER TIPS	Finnpipette F1/F2/F3 1-10µl micro	32 mm
			
KF-TFB-015L	Sterile tip with filter 15µl, 10x96 PCS/box	N/A	micro
15µl	10µl FILTER TIPS EXTENDED	Finnpipette F1/F2/F3 1-10µl micro	45 mm
			
BBA-TP-20-C-F	Sterile tip with filter 20µl, 10x96 PCS/box	Finntip Filter 20 Univ.	middle
20µl	20µl FILTER TIPS	Finnpipette F1/F2/F3 1-20µl univ	50 mm
			
BBA-BK-PT-100S	Sterile tip with filter 100µl, 10x96 PCS/box	Finntip Filter 100	middle
100µl	100µl FILTER TIPS	Finnpipette F1/F2/F3 10-100µl	50 mm
			
BK-PT-200SF	Sterile tip with filter 200µl, 10x96 PCS/box	Finntip Filter 200	middle
200µl	200µl FILTER TIPS	Finnpipette F1/F2/F3 20-200µl	51 mm
			
BK-PT-300SF	Sterile tip with filter 250µl, 10x96 PCS/box	Finntip Filter 200 Ext.	middle
250µl	200µl FILTER TIPS EXTENDED	Finnpipette F1/F2/F3 20-200µl	58 mm
			
BK-PT-1000SF	Sterile tip with filter 1000µl, 10x96 PCS/box	Finntip Filter 1000	large
1000µl	1000µl FILTER TIPS	Finnpipette F1/F2/F3 100-1000µl	78 mm
			
BK-PT-1250SF	Sterile tip with filter 1250µl, 10x96 PCS/box	Finntip Filter 1000 Ext. or Finntip Filter 1200	large
1250µl	1000µl FILTER TIPS EXTENDED	Finnpipette F1/F2/F3 100-1000µl	102 mm

VERIFIED  
BY COVID-19  
PANDEMIC:  
QUALITY &  
AVAILABILITY

